

HS-2750

Rear-Wheel Drive Foldable Powerchair

- Maximum Speed: 8 kmph / 5 mph
- Capacity: 150 kg / 331 lbs
- Lightweight aluminum frame
- Dual suspension for enhanced comfort
- Height & angle adjustable footplate
- Flick-to-fold function

TECHNICAL SPECIFICATIONS

Overall Length	1055 mm / 42" (Folded: 440 mm / 17.3")
Ovarall Width	635 mm / 25" (Folded: 635 mm / 25")
Overall Height	1040 mm / 40.9" (Folded: 890 mm / 35")
Wheels: Front	200 mm / 7.8"
Wheels: Rear	320 mm / 12.6"
Weight w/ Batteries	29.5 kgs / 65 lbs
Max. Speed	8 kmph / 5 mph
Weight Capacity	150 kg / 331 lbs
Ground Clearance	80 mm / 3.1"
Grade Climbable	6 degree
Curb Climbing	30 mm / 1.2"
Turning Radius	950 mm / 37.4"
Brake	Electro-mechanical
Seat Type	Padded Breathable Ballistic Nylon Back Rest
Seat Width	460 mm / 18.1"
Motor Size	24V / 250W / 3250 rpm
Battery Size	(1) 25.9V - 11.6Ah (25.2V - 30.15Ah Option)
Battery Weight	2.45 kgs / 5.4 lbs
Travel Range	14 km / 8.7 mi
Battery Charger	2A Off Board
Colours Available	Silver

Your local dealer is:

