



The compact Eco-Buggy is exactly right for the quick, easy transportation of children. The Eco-Buggy is especially agile and manoeuvrable.



Special features

- Higher load compared to conventional buggies
- Lightweight and very compact when folded
- Easy and uncomplicated handling



Uses

- Outdoors
- For children who are not willing or able to cover long distances on foot
- For fast, easy transportation



**ALLIED
MEDICAL**
LIMITED

Standard model

- Frame colour: anthracite metallic
- Upholstery colour: blue
- Reinforced seat and back support upholstery
- Height-adjustable, separable and swing-away leg support
- Twin-roller casters with swivel lock (pin lock)
- Wheel locks at the rear wheels
- Standard lap belt
- Folding storage bag

Options / additional accessories

- Simple padding or padding with lateral guides
- 5-point belt (only in combination with simple padding)
- Sun canopy
- Padded head support



▸ Back support angle adjustment



▸ Height-adjustable leg support, separable and swing-away



▸ Wheel lock



▸ Folding storage bag



▸ Eco-Buggy, folded



▸ Sun canopy

Technical data

Seat width	350 mm
Seat depth	300 mm
Back support height	550 mm
Back support angle	90° to 110°
Lower leg length (max. length requires caster wheel adjustment)	200–340 mm
Wheel diameter front / rear	170 mm / 210 mm
Folding size (LxVxH)	1150 x 420 x 290 mm
Weight	9 kg
Max. load	50 kg



Phone: 0800 31 61 81

Fax: 09 415 1686

Email: helpis@alliedmedical.co.nz

www.alliedmedical.co.nz