



PO Box 302 250
North Harbour
Auckland



**ALLIED
MEDICAL**
LIMITED

0800 31 61 81
www.alliedmedical.co.nz

GoTo VINYL

Our lightweight,
portable seat
for kids with
disabilities just
**got even easier
to use!**

Size 1



Promoting inclusion
and participation
in everyday life,
GoTo's new
wipe-clean vinyl
cover set make it
perfect for nursery,
school, therapy
and busy homes.

NEW FEATURES

- Wipe-clean
- Hard-wearing
- Water-resistant
- Soft-touch
- Vibrant new colour

www.fireflyfriends.com



Headrest Choice
(Standard or Advanced)

GREAT FUNCTIONS

5-point harness

Adjustable lateral support

Vinyl pressure relief cushion

Floorsitter
(can also be used with Splashy bath seat)



GoTo Size 1 and 2 suit children aged between 1-8 years (approx) who would benefit from extra postural support.

| | Seat base to top of headrest Max (mm/ inches) | Seat base to top of headrest Min (mm/ inches) | Lateral support width Max (mm/ inches) | Lateral support width Min (mm/ inches) | Seat Depth (mm/ inches) | Max Weight (kg/lbs) | Age Range (approx) |
|-----------|---|---|--|--|-------------------------------|------------------------|-----------------------|
| S1 | 540/21 | 380/14 | 260/10 | 190/7.5 | 157/7 | 15/33* | 1-3 |
| S2 | 680/27 | 510/20 | 300/12 | 210/8 | 244/9.5 | 30/66* | 3-8 |

*(if child is above 13kg/28.5lbs GoTo cannot be used in a shopping trolley)