



Scoot

Scoot is a 4-in-1 mobility rider. It lets kids with special needs play, explore and participate in ways that suit their own abilities



firefly Scooot

Scoot is a 4-in-1 mobility rider. It lets children with special needs play, explore and participate.

- ✓ Suitable for 2-6 years
- ✓ Max weight of 22kgs (44lbs)
- ✓ Can be used outside on flat, level surfaces



There are 4 configurations for Scoot



Crawl
Lying in prone, propelling
with arms and feet

Scoot
Seated and self-propelling
with feet

Ride
Seated and self-propelling
using wheels

Poosh
Seated with an adult
pushing with the handle