

Gel Ovarations™



DeZyne™ GEL

Sta-Rite™

Gel
Ovarations
Solutions



ALLIED
MEDICAL
LIMITED

UPPER BODY

Armrest Pad

GO310G (3.5" x 10")

GO312G (3.5" x 12")

GO314G (3.5" x 14")

GO316G (3.5" x 16")

GO212G (2" x 12")

- These are available in five different sizes to suit all different needs and will easily bolt onto existing armrest T sections
- They consist of a solid wood base gel pad with a neoprene nylon cover and plastic base guard
- Designed for the ultimate in comfort and protection and can aid in healing existing sores and/or abrasions
- Sold individually



GO212G

Ottobock Trough Liner

GOOTLM (Medium)

GOOTLL (Large)

- These are pads designed to fit inside Ottobock Channel Forearm Pads
- They simply Velcro in place and are long enough to pad up behind the elbow, making it considerably softer for the user
- Medium fits OB436A2=2-7 (13")
- Large fits OB436A2=1-7 (19")



LOWER BODY

Gel Knee Protector Buttons

GOKBM (Regular)
GOKBMT (Thin)

- Designed to fit swingaway legrests, these pads come complete with mounting hardware to interface with popular brands
- Regular is 1 1/4" thick, 4" diameter
- Thin is 7/8" thick, 4" diameter



Split Calf Support Panel With Strap

GOGHS2

- Split calf panel allows use of swingaway leg without removal of pad
- Fits 16" to 20" wide seats
- Positioning strap included
- Fits Right or Left leg
- Sold Individually



Leg Protector Wraps

GOLPW-LOGO

- With a velcro closure these are designed to fit around a standard 70° legrest hanger
- Sold as a pair



LOWER BODY

Footplate Protector Cover

GOGFPC

- 6" x 8" Bolt on gel footplate cover
- Designed to bolt onto standard adult sized footplates
- Wood base gel pad with neoprene cover
- Sold individually



Gel & Foam Footplate Protector

GOGFPCA1

- 6" x 8" Lightweight gel and foam footplate covers
- Designed to stick on with an adhesive back onto standard adult sized footplates
- Sold Individually



One Piece Footplate Protector Cover

GOGFPCS

- 11" x 9" gel footboard cover
- Designed to bolt onto centre-mounted footboard, common on powerchairs including the Quantum range - regular size
- It can also be used on some scooters
- Wood base gel pad with neoprene cover
- Sold individually



ADDUCTOR & LATERAL SUPPORTS

Sta-Rite Fix Bracket Adductor

GOSRABC (7/8")

GOSRABD (1")

- Quadra-pivot positioning allows pad position adjustment in 3 planes and 6 directions
- Comes standard with 3" x 4" GOAEP pad
- Serrated leg clamps and set screws at every pivot means the pad will stay right where you want it
- Size taken from 7/8" or 1" round tube leg rests
- Sold individually



Sta-Rite Swing-Away Bracket Adductor

GOSRASC (7/8")

GOSRASD (1")

GOSRASF (Flat)

- Quadra-pivot positioning allows pad position adjustment in 3 planes and 6 directions
- Comes standard with 3" x 4" GOAEP pad
- Serrated leg clamps and set screws at every pivot means the pad will stay right where you want it
- 7/8" or 1" fits round tube leg rests
- FM fits flat to square tube seating systems or surfaces
- Sold as a pair



Sta-Rite 6-Way Lateral Thigh Support

GOSRMHG10 (3.5" x 10")

GOSRMHG12 (3.5" x 12")

GOSRMHG14 (3.5" x 14")

- Multi-axis Hip Guide that is adjustable in 6 directions
- They mount easily to the side of virtually any seating system and can be easily adjusted up and down and can be locked in the position of choice
- Made of stainless steel on the primary support part with 2 interlocking crown pivots to allow adjustment in and out. The track on the pad allows continuous forward and back adjustment
- Sold individually



PADS

Oval Stick On Pad

GO34SP (3" x 4")

GO48SP (4" x 8")

GO412SP (4.5" x 12")

- Flat stick on gel pad with neoprene nylon cover
- Each is 0.375" in thickness
- Ideal to place on trays, calf pads and any surface that the user's body may come into contact with
- Sold individually



Egg Shaped Stick On Pad

GOESPLOGO (2.625" x 4.5")

- Domed stick on gel pad with neoprene nylon cover
- Each is 0.625" in thickness
- Sold individually



Lateral Gel Pad

GO46DGLP1

- 4" x 6" Wheelchair Lateral Pad made with silicone Dimensional gel
- Removable washable Neoprene / Nylon cover
- Standard 1" T-nut spacing fits most lateral support hardware
- Sold individually



Adductor / Elbow Gel Pad

GOAEP


- Flat 3" x 4" gel pad comes with a solid wooden bottom
- Bottom has a plastic cover and 2 T-nut connections with 1" spacing
- Sold individually



DEZYNE GEL SHEETS & PADS

DeZyne Gel Sheet

GODG14145 (14" x 14")
GODG16165 (16" x 16")
GODG16185 (16" x 18")
GODG18165 (18" x 16")
GODG18185 (18" x 18")
GODG18205 (18" x 20")
GODG20185 (20" x 18")
GODG20205 (20" x 20")
GODG22185 (22" x 18")
GODG22205 (22" x 20")
GODG24205 (24" x 20")



- The unique lightweight cellular design can redistribute pressure, reduce abrasion, absorb light shocks and dampen vibration while minimising heat and moisture build up
- Made from water insoluble proprietary silicone gel so that it is ideal for use in places where the pad might get wet, so it's easy to clean
- Adapts to all surfaces and is long lasting as it will not flow or leak if cut or punctured, and is also odourless, hypoallergenic and latex free
- Each is 0.5" in thickness

DeZyne Gel Pad

GODG355 (3" x 5")
GODG445 (4" x 4")
GODG45755 (4.5" x 7.5")
GODG465 (4" x 6")
GODG555 (5" x 5")
GODG57555 (5" x 7")
GODG65955 (6.5" x 9.5")
GODG775 (7" x 7")
GODG895 (8" x 9")



- Pad sized version of the DeZyne Gel Sheet
- The unique lightweight cellular design can redistribute pressure, reduce abrasion, absorb light shocks and dampen vibration while minimising heat and moisture build up
- Made from water insoluble proprietary silicone gel so that it is ideal for use in places where the pad might get wet, so it's easy to clean
- Long lasting, will not flow or leak if cut or punctured, it is also odourless, hypoallergenic and latex free
- Adapts to all surfaces
- Each is 0.5" in thickness

ACCESSORIES

DeZyne Gel Sheet Cover

ROLAGELCOVxxxx

- Available in sizes to match the DeZyne Gel Sheets (14" x 14" up to 24" x 20")
- This is to keep the gel sheets clean and easy to care for
- Made of black urethane fabric that is easy to wipe down
- Has a zip at the back for easy removal
- Sold individually



Screw Sets for Gel Ovation Armpads

GORTFSL

- Adapter set used to convert 3/4" round tube to flat surface using 2 bolts
- Most commonly used with flat bottom arm pads
- The 2 adapter pieces accommodate various bolt configurations
- Order 2 sets to adapt a pair of flat bottom arm rests



Serrated Edge Tube Clamp

GOTC1 (1")

GOTC2 (2")

- Both sizes of the serrated edge tube clamp will fit 7/8" and 1" diameter tubing but not 3/4"
- Sold individually



www.alliedmedical.co.nz
PO Box 302 250, North Harbour, Auckland
T: 0800 31 61 81 E: sales@alliedmedical.co.nz

**ALLIED
MEDICAL**
LIMITED