



HS-2850

Compact Mid-Wheel Drive Powerchair

- Maximum Speed : 6.7 kmph / 4.2 mph
- Capacity: 136 kg / 300 lbs
- Strong and responsive front and rear castors
- Length adjustable joystick
- Optional swing-away joystick
- Optional elevating legrests

TECHNICAL SPECIFICATIONS

Overall Length	995 mm / 39"
Overall Width	610 mm / 24"
Overall Height	1080 mm / 42.5"
Wheels: Front	155 mm / 6"
Wheels: Middle	260 mm / 10"
Wheels: Rear	155 mm / 6"
Weight w/ Batteries	82 kg / 181 lbs
Max. Speed	6.7 kmph / 4.2 mph
Weight Capacity	136 kg / 300 lbs
Ground Clearance	60 mm / 2.4"
Grade Climbable	6 degree
Curb Climbing	45 mm / 1.8"
Turning Radius	620 mm / 24.4"
Brake	Electro-mechanical
Seat Type	Swivel mid-back w/ headrest, adjustable seat height & seat belt
Seat Width	490 mm / 19"
Motor Size	420W, 4900 rpm
Battery Size	(2) 12V. 36Ah
Battery Weight	22.5 kg / 49.6 lbs
Travel Range	20 km / 12.4 miles (dependent on user weight and terrain)
Battery Charger	5A Off board
Colours Available	Red



Your local dealer is: