



HS-328

Mid-Range Scooter

- Maximum Speed : 8 kmph / 5 mph
- Capacity : 136 kg / 300 lbs
- Easy Breakdown Mechanism
- Patented Battery Pack Locking Mechanism
- Standard LED Headlight
- 22Ah Battery Pack For Longer Travel Range
- Standard Delta Handle Bar
- Optional Tiller Head Charging Socket

TECHNICAL SPECIFICATIONS

Overall Length	1120 mm / 44"
Overall Width	590 mm / 23"
Overall Height	935 mm / 37"
Wheels: Front	225 mm / 9"
Wheels: Rear	225 mm / 9"
Weight w/ Batteries	55 kg / 121 lbs
Weight of Heaviest Piece	17.4 kg / 38lbs
Max. Speed	8 kmph / 5 mph
Weight Capacity	136 kg / 300 lbs
Ground Clearance	52 mm / 2"
Grade Climbable	8 degree
Curb Climbing	50 mm / 2"
Turning Radius	1450 mm / 57"
Suspension	Rear (Option)
Brake	Electro-Mechanical
Seat Type	Swivel Light Weight Foldable
Seat Width	460 mm / 18"
Motor Size	250W, 4000 r.p.m
Battery Size	(2) 12V. 22Ah
Battery Weight	15.5 kg / 34 lbs
Travel Range	17.5 km / 11 Miles (dependent on user weight and terrain)
Battery Charger	3A Off Board, Battery Level Indicator,
Electronics	On / Off Key Switch, Speed Control Knob
Colours Available	Blue

Your local dealer is:

