



HS-559

Delux Mid-Range Scooter

- Maximum Speed : 12 kmph / 7.5 mph
- Capacity : 136 kg / 300 lbs
- 290 mm / 11" Front And Rear Wheels
- Premium Upholstery Mid-back Swivel Seat
- Seat Sliding Adjustment
- Integrated Front and Rear Bumper
- Anti Fatigue Delta Handle Bar
- Full Lighting Package
- Full Coil-Over Suspension

TECHNICAL SPECIFICATIONS

Overall Length	1200 mm / 47.2"
Overall Width	628 mm / 24.7"
Overall Height	1190 mm / 47"
Wheels: Front	290 mm / 11"
Wheels: Rear	290 mm / 11"
Weight w/ Batteries	101 kg / 222 lbs
Max. Speed	12 kmph / 7.5 mph
Weight Capacity	136 kg / 300 lbs
Ground Clearance	80 mm / 3.1"
Grade Climbable	10 degree
Curb Climbing	50 mm / 2"
Turning Radius	1500 mm / 59.1"
Suspension	Front & Rear
Brake	Electro-Mechanical
Seat Type	Premium Upholstery Mid-Back Swivel Sliding Seat
Seat Width	455 mm / 18"
Motor Size	700W 4500 rpm
Battery Size	(2) 12V. 50Ah
Battery Weight	30.5 kg / 67.2 lbs
Travel Range	42.5 km / 26.4 miles (dependent on user weight and terrain)
Battery Charger	5A Off Board, Battery Level Indicator,
Electronics	On /Off Key Switch, Speed Control Knob
Colours Available	Blue, Burgundy, Silver



Your local dealer is: