



HS-5600

Heavy Duty Rear-Wheel Drive Powerchair

- Maximum Speed : 8 kmph / 5 mph
- Weight capacity: 181 kg / 400 lbs
- Effortless single free-wheel lever
- Vertical and horizontal adjustable footplate
- Swing-away joystick
- Optional elevating legrests

TECHNICAL SPECIFICATIONS

Overall Length	1135 mm / 44.7"
Overall Width	635 mm / 25"
Overall Height	1185 mm / 46.7"
Wheels: Front	230 mm / 9"
Wheels: Rear	360 mm / 14"
Weight w/ Batteries	102 kg / 224 lbs
Max. Speed	8 kmph / 5 mph
Weight Capacity	181 kg / 400 lbs
Ground Clearance	85 mm / 3.4"
Grade Climbable	8 degree
Curb Climbing	60 mm / 2.4"
Turning Radius	790 mm / 31.1"
Brake	Electro-mechanical
Seat Type	Captain seat, adjustable seat height and seat belt
Seat Width	500 mm / 20"
Motor Size	420W, 4600 rpm
Battery Size	(2) 12V. 50Ah
Battery Weight	30.5 kg / 67.2 lbs
Travel Range	24 km / 15 miles (dependent on user weight and terrain)
Battery Charger	5A Off board
Colours Available	Red



Your local dealer is: