



# HS-898

## Bravo Edition

### Hill Climber Scooter

- Maximum Speed : 10 kmph / 6.2 mph
- Capacity : 181 kg / 399 lbs
- 330 mm / 13" Front and Rear Wheels
- Anti Fatigue Delta Handle Bar
- Full Lighting Package
- Full Coil-Over Suspension
- Hand Brake
- Hi / Lo Speed Switch

## TECHNICAL SPECIFICATIONS

Overall Length	1440 mm / 57"
Overall Width	720 mm / 28.4"
Overall Height	1360 mm / 53.5"
Wheels: Front	330 mm / 13"
Wheels: Rear	330 mm / 13"
Weight w/ Batteries	143.4 kg / 316.1 lbs
Max. Speed	10 kmph / 6.2 mph
Weight Capacity	181 kg / 400 lbs
Ground Clearance	100 mm / 3.9"
Grade Climbable	10 degree
Curb Climbing	80 mm / 3.2"
Turning Radius	1600 mm / 63"
Suspension	Front & Rear
Brake	Manual Hand Brake & Electro-Mechanical
Seat Type	Sliding and Swivel Reclining Captain W/Headrest
Seat Width	520 mm / 20.5"
Motor Size	850W, 3450 rpm
Battery Size	(2) 12V. 75Ah
Battery Weight	55 kg / 121 lbs
Travel Range	40 km / 25 miles (dependent on user weight and terrain)
Battery Charger	8A Off Board
Electronics	On /Off Key Switch, Speed Control Knob
Colours Available	Blue, Burgundy

Your local dealer is:

