



# HS-898

## Heavy Duty Scooter

- Maximum Speed : 15 kmph / 9.3 mph
- Capacity : 181 kg / 400 lbs
- 330 mm / 13" Front and Rear Wheels
- Anti Fatigue Delta Handle Bar
- Full Lighting Package
- Reclining Captains Chair
- Full Coil-Over Suspension
- Hand Brake
- Hi / Lo Speed Switch

### TECHNICAL SPECIFICATIONS

Overall Length	1440 mm / 57"
Overall Width	680 mm / 26.8"
Overall Height	1360 mm / 53.5"
Wheels: Front	330 mm / 13"
Wheels: Rear	330 mm / 13"
Weight w/ Batteries	117 kg / 258 lbs
Max. Speed	15 kmph / 9.3 mph
Weight Capacity	181 kg / 400 lbs
Ground Clearance	100 mm / 3.9"
Grade Climbable	10 degree
Curb Climbing	80 mm / 2.8"
Turning Radius	1600 mm / 63"
Suspension	Front & Rear
Brake	Manual Hand Brake & Electro-Mechanical
Seat Type	Sliding and Swivel Reclining Captain W/Headrest
Seat Width	520 mm / 20.5"
Motor Size	850W, 5000 rpm
Battery Size	(2) 12V. 50Ah
Battery Weight	30.5 kg / 67.2 lbs
Travel Range	27 km / 16.8 miles (dependent on user weight and terrain)
Battery Charger	5A Off Board, Battery Level Indicator
Electronics	On /Off Key Switch, Speed Control Knob
Colours Available	Burgundy



Your local dealer is: