



HS-928

Deluxe Long Range Scooter with Cabin

- Maximum Speed : 10 kmph / 6.2 mph
- Capacity : 200 kg / 441 lbs
- Front windscreen wiper
- Rearview mirror
- Tail lights and rear turning signal
- Two 12V.100Ah Batteries for longer travel range
- Swivel/Sliding Captains Chair w/ angle adjust
- Heater and indoor light

TECHNICAL SPECIFICATIONS

Overall Length	1600 mm / 63"
Overall Width	1040 mm / 40.9"
Overall Height	1620 mm / 63.8"
Wheels: Front	380 mm / 15"
Wheels: Rear	380 mm / 15"
Weight w/ Batteries	229 kg / 505 lbs
Max. Speed	10 kmph / 6.2 mph
Weight Capacity	200 kg / 441 lbs
Ground Clearance	110 mm / 4.3"
Grade Climbable	10 degree
Curb Climbing	100 mm / 3.9"
Turning Radius	1860 mm / 73.2"
Suspension	Front & Rear
Brake	Manual Hand Brake & Electro-Mechanical
Seat Type	Sliding and Swivel Captain Style w Back Angle Adjustment
Seat Width	559 mm / 22"
Motor Size	1000W, 3200 rpm
Battery Size	(2) 12V. 100Ah
Battery Weight	65.3 kg / 144 lbs
Travel Range	44 km / 27 miles (dependent on user weight and terrain)
Battery Charger	12A Off Board, Battery Level Indicator
Electronics	On /Off Key Switch, Speed Control Button
Colours Available	Silver



Your local dealer is: