



HS-928

Long Range Scooter

- Maximum Speed : 15 kmph / 9.3 mph
- Capacity : 226 kg / 500 lbs
- LCD Display
- Tiller Storage Compartment
- Large 15" Front & Rear Wheels
- Two 12V/100Ah Batteries for longer travel range
- Reclining Captains Chair w/ angle adjust seat

TECHNICAL SPECIFICATIONS

Overall Length	1600 mm / 63"
Overall Width	720 mm / 28.3"
Overall Height	1280 mm / 50"
Wheels: Front	380 mm / 15"
Wheels: Rear	380 mm / 15"
Weight w/ Batteries	178 kg / 392 lbs
Max. Speed	15 kmph / 9.3 mph
Weight Capacity	226 kg / 500 lbs
Ground Clearance	110 mm / 4.3"
Grade Climbable	9 degree
Curb Climbing	90 mm / 3.5"
Turning Radius	1860 mm / 73.2"
Suspension	Front & Rear
Brake	Manual Hand Brake & Electro-Mechanical
Seat Type	Sliding and Swivel Reclining Captain W/Headrest
Seat Width	559 mm / 22"
Motor Size	800W, 4300 rpm
Battery Size	(2) 12V. 100Ah
Battery Weight	65.3 kg / 144 lbs
Travel Range	50 km / 31 miles (dependent on user weight and terrain)
Battery Charger	12A Off Board, Battery Level Indicator
Electronics	On /Off Key Switch, Speed Control Button
Colours Available	Blue, Burgundy

Your local dealer is:

