




# Ride® Java® Back Regular Size Matrix



**REGULAR**



**DEEP**

Item	Measurement Type
<b>A</b>	Overall Back Width
<b>B</b>	Width Inside Foam
<b>C</b>	Back Height
<b>D</b>	Back Depth





<b>REGULAR</b>		<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>		
Back Width	Fits Wheelchair Widths	Overall Back Width	Width (Inside Foam)	Back Height	Back Depth	Carrying Weight*	On-Chair Weight**
14" Regular	13-16" / 33-41cm	14" / 36cm	13" / 33cm	14" / 36cm	3" / 8cm	2.0 lb / 1.0kg	3.8 lb / 1.7kg
15" Regular	14-17" / 36-43cm	15" / 38cm	13.5" / 34cm	14" / 36cm	3" / 8cm	2.4 lb / 1.1kg	4.2 lb / 1.9kg
16" Regular	15-18" / 38-46cm	16" / 41cm	14.5" / 37cm	14" / 36cm	3" / 8cm	2.7 lb / 1.2kg	4.5 lb / 2.0kg
17" Regular	16-19" / 41-48cm	17" / 43cm	15.5" / 39cm	14" / 36cm	3" / 8cm	2.9 lb / 1.3kg	4.7 lb / 2.1kg
18" Regular	17-20" / 43-51cm	18" / 46cm	16.5" / 42cm	14" / 36cm	3" / 8cm	3.1 lb / 1.4kg	4.9 lb / 2.2kg
19" Regular	18-21" / 46-53cm	19" / 48cm	17.5" / 44cm	14" / 36cm	3" / 8cm	3.3 lb / 1.5kg	5.1 lb / 2.3kg
20" Regular	19-22" / 48-56cm	20" / 51cm	18.5" / 47cm	14" / 36cm	3" / 8cm	3.5 lb / 1.6kg	5.3 lb / 2.4kg

<b>REGULAR DEEP</b>		<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>		
Back Width	Fits Wheelchair Widths	Overall Back Width	Width (Inside Foam)	Back Height	Back Depth	Carrying Weight*	On-Chair Weight**
14" Regular	13-16" / 33-41cm	14" / 36cm	13" / 33cm	14" / 36cm	6" / 15cm	2.4 lb / 1.1kg	4.2 lb / 1.9kg
15" Regular	14-17" / 36-43cm	15" / 38cm	13.5" / 34cm	14" / 36cm	6" / 15cm	2.7 lb / 1.2kg	4.5 lb / 2.0kg
16" Regular	15-18" / 38-46cm	16" / 41cm	14.5" / 37cm	14" / 36cm	6" / 15cm	3.0 lb / 1.3kg	4.8 lb / 2.2kg
17" Regular	16-19" / 41-48cm	17" / 43cm	15.5" / 39cm	14" / 36cm	6" / 15cm	3.2 lb / 1.4kg	5.0 lb / 2.2kg
18" Regular	17-20" / 43-51cm	18" / 46cm	16.5" / 42cm	14" / 36cm	6" / 15cm	3.4 lb / 1.5kg	5.2 lb / 2.3kg
19" Regular	18-21" / 46-53cm	19" / 48cm	17.5" / 44cm	14" / 36cm	6" / 15cm	3.6 lb / 1.6kg	5.4 lb / 2.4kg
20" Regular	19-22" / 48-56cm	20" / 51cm	18.5" / 47cm	14" / 36cm	6" / 15cm	3.8 lb / 1.7kg	5.6 lb / 2.5kg

Specifications are subject to change. \* Java Back, less mounting hardware, \*\* Java Back, complete with hardware.



# Ride® Java® Back Tall Size Matrix



Item	Measurement Type
<b>A</b>	Overall Back Width
<b>B</b>	Width Inside Foam
<b>C</b>	Back Height
<b>D</b>	Back Depth



<b>TALL</b>		<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>		
Back Width	Fits Wheelchair Widths	Overall Back Width	Width (Inside Foam)	Back Height	Back Depth	Carrying Weight*	On-Chair Weight**
14" Tall	13-16" / 33-41cm	14" / 36cm	13" / 33cm	18" / 46cm	3" / 8cm	2.7 lb / 1.2kg	4.5 lb / 2.1kg
15" Tall	14-17" / 36-43cm	15" / 38cm	13.5" / 34cm	18" / 46cm	3" / 8cm	2.9 lb / 1.3kg	4.7 lb / 2.2kg
16" Tall	15-18" / 38-46cm	16" / 41cm	14.5" / 37cm	18" / 46cm	3" / 8cm	3.2 lb / 1.5kg	5.0 lb / 2.3kg
17" Tall	16-19" / 41-48cm	17" / 43cm	15.5" / 39cm	18" / 46cm	3" / 8cm	3.4 lb / 1.6kg	5.2 lb / 2.4kg
18" Tall	17-20" / 43-51cm	18" / 46cm	16.5" / 42cm	18" / 46cm	3" / 8cm	3.6 lb / 1.7kg	5.4 lb / 2.5kg
19" Tall	18-21" / 46-53cm	19" / 48cm	17.5" / 44cm	18" / 46cm	3" / 8cm	3.8 lb / 1.7kg	5.6 lb / 2.6kg
20" Tall	19-22" / 48-56cm	20" / 51cm	18.5" / 47cm	18" / 46cm	3" / 8cm	4.0 lb / 1.8kg	5.8 lb / 2.7kg

<b>TALL DEEP</b>		<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>		
Back Width	Fits Wheelchair Widths	Overall Back Width	Width (Inside Foam)	Back Height	Back Depth	Carrying Weight*	On-Chair Weight**
14" Tall	13-16" / 33-41cm	14" / 36cm	13" / 33cm	18" / 46cm	6" / 15cm	3.3 lb / 1.5kg	5.1 lb / 2.3kg
15" Tall	14-17" / 36-43cm	15" / 38cm	13.5" / 34cm	18" / 46cm	6" / 15cm	3.3 lb / 1.5kg	5.1 lb / 2.3kg
16" Tall	15-18" / 38-46cm	16" / 41cm	14.5" / 37cm	18" / 46cm	6" / 15cm	3.6 lb / 1.7kg	5.4 lb / 2.5kg
17" Tall	16-19" / 41-48cm	17" / 43cm	15.5" / 39cm	18" / 46cm	6" / 15cm	3.8 lb / 1.7kg	5.6 lb / 2.6kg
18" Tall	17-20" / 43-51cm	18" / 46cm	16.5" / 42cm	18" / 46cm	6" / 15cm	4.0 lb / 1.8kg	5.8 lb / 2.7kg
19" Tall	18-21" / 46-53cm	19" / 48cm	17.5" / 44cm	18" / 46cm	6" / 15cm	4.2 lb / 2.0kg	6.0 lb / 2.7kg
20" Tall	19-22" / 48-56cm	20" / 51cm	18.5" / 47cm	18" / 46cm	6" / 15cm	4.4 lb / 2.0kg	6.2 lb / 2.8kg

Specifications are subject to change. \* Java Back, less mounting hardware, \*\* Java Back, complete with hardware.©



# Ride® Java® Back Grande Size Matrix



Item	Measurement Type
<b>A</b>	Overall Back Width
<b>B</b>	Width Inside Foam
<b>C</b>	Back Height
<b>D</b>	Back Depth

GRANDE REGULAR		A	B	C	D	Carrying Weight*	On-Chair Weight**
Back Width	Fits Wheelchair Widths	Overall Back Width	Width (Inside Foam)	Back Height	Back Depth		
14" Grande	13-16" / 33-41cm	14" / 36cm	13" / 33cm	22" / 56cm	3" / 8cm	3.0 lb / 1.4kg	4.8 lb / 2.2kg
15" Grande	14-17" / 36-43cm	15" / 38cm	13.5" / 34cm	22" / 56cm	3" / 8cm	4.0 lb / 1.8kg	5.8 lb / 2.7kg
16" Grande	15-18" / 38-46cm	16" / 41cm	14.5" / 37cm	22" / 56cm	3" / 8cm	4.4 lb / 2.0kg	6.2 lb / 2.8kg
17" Grande	16-19" / 41-48cm	17" / 43cm	15.5" / 39cm	22" / 56cm	3" / 8cm	4.6 lb / 2.1kg	6.4 lb / 2.9kg
18" Grande	17-20" / 43-51cm	18" / 46cm	16.5" / 42cm	22" / 56cm	3" / 8cm	4.8 lb / 2.2kg	6.6 lb / 3.0kg
19" Grande	18-21" / 46-53cm	19" / 48cm	17.5" / 44cm	22" / 56cm	3" / 8cm	5.0 lb / 2.3kg	6.8 lb / 3.1kg
20" Grande	19-22" / 48-56cm	20" / 51cm	18.5" / 47cm	22" / 56cm	3" / 8cm	5.2 lb / 2.4kg	7.0 lb / 3.2kg

GRANDE DEEP		A	B	C	D	Carrying Weight*	On-Chair Weight**
Back Width	Fits Wheelchair Widths	Overall Back Width	Width (Inside Foam)	Back Height	Back Depth		
14" Grande	13-16" / 33-41cm	14" / 36cm	13" / 33cm	22" / 56cm	6" / 15cm	3.3 lb / 1.5kg	5.1 lb / 2.3kg
15" Grande	14-17" / 36-43cm	15" / 38cm	13.5" / 34cm	22" / 56cm	6" / 15cm	4.4 lb / 2.0kg	5.2 lb / 2.4kg
16" Grande	15-18" / 38-46cm	16" / 41cm	14.5" / 37cm	22" / 56cm	6" / 15cm	4.8 lb / 2.2kg	6.6 lb / 3.0kg
17" Grande	16-19" / 41-48cm	17" / 43cm	15.5" / 39cm	22" / 56cm	6" / 15cm	5.0 lb / 2.3kg	6.8 lb / 3.1kg
18" Grande	17-20" / 43-51cm	18" / 46cm	16.5" / 42cm	22" / 56cm	6" / 15cm	5.2 lb / 2.4kg	7.0 lb / 3.2kg
19" Grande	18-21" / 46-53cm	19" / 48cm	17.5" / 44cm	22" / 56cm	6" / 15cm	5.4 lb / 2.5kg	7.2 lb / 3.3kg
20" Grande	19-22" / 48-56cm	20" / 51cm	18.5" / 47cm	22" / 56cm	6" / 15cm	5.6 lb / 2.6kg	7.4 lb / 3.4kg

Specifications are subject to change. \* Java Back, less mounting hardware, \*\* Java Back, complete with hardware.