

Client: _____

Assessor: _____

Organisation: _____

Karma BT10 Bariatric Wheelchair Script

Please select one option from each available or suitable section



STANDARD FEATURES

- Tension adjustable back upholstery
 - Solid seat base
 - Swingaway 70 degree legrests with angle adjustable footplates
 - Flip back armrests - fixed height
 - Angle adjustable back upholstery (from -5° to +15°)
 - 2 point hip belt
 - Heel loops
 - Swivel anti-tippers
 - Basic foam cushion
 - Dual brake system - attendant brake as well as manual brake on each wheel
 - Rear wheel 24" pneumatic
 - Castor 7" x 1 3/4"
 - Standard seat to floor height 17" (is adjustable)
- ‡ = Standard Feature
 Ø = Free of Charge
 \$ = Extra Cost Option

SEAT WIDTH AND SEAT DEPTH	WIDTH			
	558 mm (22")	609 mm (24")	660 mm (26")	
457 mm (18")				Ø
508 mm (20")				Ø

Armrest Type

Seat to Floor Height

Rear Wheels

Rear Axle

Space between Pushrim
& Rear Wheel

Seat Positioning Belt

Back Accessories

Accessories

IV Pole (KMIVPOLE)

Oxygen Bottle Holder (KMO2)

Tray - slide on (TRAYKMACRXSO)

Legrest Type

Front Castor

Brakes

Backrest Height (adjustable)

Back Angle (adjustable)

Seat Accessories

Notes