

Karma Flexx Wheelchair Script



Customer Details							
Client							
Assessor							
Organisation							
	ACC			MOH		Other	
User Weight							

Karma Flexx Wheelchair

Please note: Scripts are not required for Karma Flexx wheelchairs with the standard specifications (listed below). This script should only be used for clients with alternative requirements.

Standard Specifications:

- 7" Castor Wheels
- 19" Seat to Floor Height
- 1" Seat Dump
- Height Adjustable Ergo Flipback Arms
- 80° Swing-away Leg supports
- Flip-up Foot Supports
- Push to Lock Wheel Locks
- Swing Away Anti-Tippers
- Calf Strap
- Standard Hip Belt
- Seat widths 14-18": 130kg max user weight
- Seat widths 20"+: 170kg max user weight

Wheelchair Type

- Self Propelling
- Quick Release 24" Pneumatic Spoked Rear Wheel
 - Ergo Handrim

- Transit
- Fixed 14" Solid Tyre Plastic Rim Rear Wheels

Seat Size (Width x Depth)

14 x 16
16 x 16
16 x 18
16 x 20
18 x 16
18 x 18
18 x 20

20 x 18 [^]
20 x 20 [^]
22 x 18 [^]
22 x 20 [^]
24 x 18 [^]
24 x 20 [^]

[^] Heavy Duty Model

Front Seat Height

		Seat to Floor Height					
Castor Size	1) 5"		16"	17"	18"		
	6"			17"	18"	19"	
	7"			17.5"	18"	19"	
	8"					19"	20"

¹⁾ 5" Castor not available on models 20" and wider

Back Cane Height

		Back Height					
Fixed Height	Std			17"			*
Height Adjustable (Adjustable with tools only)	\$	15"	16"	17"	18"	19"	20"

* Standard on 16 x 20 and 18 x 20 models

Rear Seat Height

		Rear Seat to Floor Height					
Rear Wheel Size	14"					18"	
	20"	14"	15"	16"	17"	18"	
	22"		15"	16"	17"	18"	19"
	24"			16"	17"	18"	19"

Anti Tippers

Swing Away Height Adjustable

Front Seat to Floor Height 18-21".....Std
Front Seat to Floor Height 15-17".....\$

- Pneumatic tyres Solid tyres

Centre of Gravity

- Rear Centre Forward

Seat Upholstery

Standard (Padded)	Std
Standard with velcro*	\$
Omit Seat Upholstery	N/C
Cushion (please specify below)	\$
Spare Cushion Cover	\$
Cushion	

*Available on most sizes. Some sizes may be indent only.

Transit Option

Transit Tie Down Option	\$
-------------------------------	----

* Not available on HD

Leg Support

	Left	Right	
80° Swing-away			Std
90° Swing-away			\$
Elevating			\$
Stump Support			\$

Foot Support

	Left	Right	
Standard Flip-up			Std
Flip-up Angle Adjustable*			\$
Flip-up Angle & Depth Adj. Small..			\$
Large (16" or wider chairs only)....			\$

* Not available on HD

Calf Strap/Heel Loops

Calf Strap	Std
Heel Loops*	\$

*Standard on Standard Flip-up only, otherwise an upcharge.

Other Parts Required

Back Support

Standard (Padded)	Std
Tension Adjustable (not available on height adj back canes)	\$
Omit Back Upholstery	N/C
Alternative Back Support (please specify below)	\$
Alternative Back Support	

Arm Support

Flip-back Ergo (desk length) Height Adjustable	Std
Flip-back Boat (full length) Height Adjustable	\$

Wheel Locks

Push to Lock with Brake extension	Std
Pull to Lock*	\$

* Available only on 20, 22 and 24" Rear wheels

Accessories

Clear Spoke Guards	\$
Oxygen Tank Holder	\$
Push Handle Rigidiser Bar*	\$
IV Pole	\$

* Std on height adj back canes 18" and wider and HD chairs

Head Support

Foldable Medium (16-18" widths)	\$
Foldable Large (20" widths)	\$
Other	\$

Pelvic Positioning Belt

Lap Belt	Std
2 Point	\$
4 Point	\$

Tech Notes