

◆ *Stand up again, improve your independence
and be kind to your body!*

KARMA
home & health care

LIVE
YOUR
LIFE


< ERGO STAND >




S-Ergo

😊 S-ERGO SEAT

PATENTED
PATENTED
PATENTED



☹️ TRADITIONAL SEAT



< KP 80 >

◆ The **Ergo Stand** is a highly manoeuvrable standing power wheelchair with *Karma's* unique, patented *S-Ergo* seating system. The front wheel drive *Ergo Stand* is only 62 cm wide and has excellent indoor manoeuvrability performance. Agile control and pleasant driving characteristics during outdoor use is guaranteed. The *Ergo Stand* integrates many adjustable functions to meet individual user requirement. The power recline function provides a wide angle adjustment of 80° – 120° with excellent back support.

- ◆ *Power stand*
- ◆ *Power backrest recline*



Features



• FRONT WHEEL DRIVE



• POWER SIT TO STAND



• POWER RECLINE



• FLIP BACK AND HEIGHT ADJUSTABLE ARMREST



• LIGHTS



• KNEE SUPPORT



• SWING AWAY JOYSTICK



Reduced sliding



Pressure redistribution

• S-ERGO SEAT SYSTEM

< The medical benefits of standing >

- ◆ Activates circulation system
- ◆ Improves cardio-pulmonary function
- ◆ Helps prevent decubitus sores
- ◆ Improves bowel regularity and function
- ◆ Helps bladder and intestinal motility
- ◆ Reduces contractures and ankylosis

SPECIFICATION

STANDARD

Speed: 6-9 km/h
 Width: 63 cm
 Length: 102 cm
 Seat Width: 40 / 44 / 48 cm
 Motor: 2x 300W
 Battery: 2x 36 Ah / 2x 50 Ah
 Range: 28 km (50 Ah)
 Maximum User Weight: 135 kg
 Wheelchair Weight: 89-91 kg
 Weight Without Batteries: 57 kg

S-Ergo seat system

Power stand

Power recline backrest

Flip back and height adjustable armrests

Knee support adjustable in height, depth and width

Light kit

Front wheel drive

Color: Diamond blue

Driving in standing position

± 1cm ± 0,5 kg

OPTIONS

Parallel swing-away joystick	Seat belt
Head support	Mud guard
H-shaped seat belt	

PO Box 302 250
 North Harbour
 Auckland

0800 31 61 81
www.alliedmedical.co.nz



**ALLIED
 MEDICAL**
 LIMITED