

Karma Wheelchair Quick Reference Guide



KM8020



KM8520



Flexx Junior



Flexx



Flexx Tall



Flexx HD

User Weight Limits

Model	Seat size	User Weight Limit
KM8020/8520	14" – 18" width	130 kg
KM8020/8520	20" – 22" width	160 kg
Flexx	14" – 18" width	130 kg
Flexx Tall	16" - 18" width	130 kg
Flexx HD	20" - 24" width	170 kg
Flexx Junior	10" – 14" width	75 kg

Flexx Range Seat to Floor Height

Rear wheel	5" castor	6" castor	7" castor	8" castor
14" (Transit)	18"	18"	18"	N/A
20"	16" - 17"	17"	N/A	N/A
22"	16" - 18"	17" - 18"	18"	N/A
24"	17" - 18"	17" - 19"	18" - 19"	19" - 21"

KM8020/8520 Seat to Floor Height**

Rear wheel	6" castor	7" castor	8" castor
12" (Transit)	N/A	N/A	19"
20"	16" * or 17"	18"	17"
22"	18"	19"	18"
24"	19"	19" / 20"	19" / 20" / 21"



0800 31 61 81
www.alliedmedical.co.nz

***With Drop Seat Base
**Seat to Floor Height measurement is taken from the front of the seat upholstery
(These measurements are a guide only and do not factor in specifics such as seat rake)*

Backrest Heights

Model	Standard Height (seat to top of back upholstery)	Adjustment available
KM8020/8520	15.5" (393 mm)	Up to 18.5" (470 mm) 22" (559 mm) with tall backcane option
Flexx / Tall	17" (432 mm)	15" - 21" (381 - 533 mm)
Flexx HD	17" (432 mm)	15" - 20" (381 - 508 mm)
Flexx Junior	14.5" (370 mm)	14.5" - 17" (370 - 430 mm)

Armrest Heights

Model	Minimum Height	Maximum Height
KM8020/8520	11.5" (292 mm)	18.5" (470 mm)
Flexx / Tall	8" (203 mm)	11.5" (292 mm)
Flexx HD	8" (203 mm)	11.5" (292 mm)
Flexx Junior	6" (150 mm)	10" (250 mm)

Armrest Lengths

Model	Desk	Full length
KM8020/8520	10" (254 mm) x 1 7/8" (47 mm)	12" (304 mm) x 1 7/8" (47 mm)
Flexx / Tall	One size - 11.5" (292 mm) x 2" (50 mm)	
Flexx HD	One size - 11.5" (292 mm) x 2" (50 mm)	
Flexx Junior	One size - 9.6" (244.5 mm) x 1.76 (44.7 mm)	

Lower Leg Length

Model	Standard Length Range	Shortest Length Achievable	Longest Length
KM8020/8520	16" - 18"(406 mm - 457 mm)	13" (330 mm)	20" (508 mm)
Flexx / Tall	15.5" - 19" (393 mm - 483mm)	15" (381 mm)	19.5" (495 mm)
Flexx HD	15.5" - 19" (393 mm - 483mm)	15" (381 mm)	19.5" (495 mm)
Flexx Junior	7" - 11" (180mm - 280mm)	7" (180 mm)	16" (400 mm)

Tube Diameter

Armrest Tube Size	3/4"
Back Cane Tube Size	7/8"

Overall Width

Model	Total width
KM8020 (Self propelling)	Seat width plus 9" (229 mm)
KM8520 (Transit)	Seat width plus 6" (152 mm)
Flexx	Seat width plus 9" (229 mm)
Flexx Tall	Seat width plus 9" (229 mm)
Flexx HD	Seat width plus 9" (229 mm)
Flexx Junior	Seat width plus 7.5" or 7.9" (190 or 200mm)

This information has been prepared to make prescribing Karma wheelchairs easier. If there are any other specific questions you have, please contact your Business Territory Manager or our Customer Care Team on **0800 31 61 81**.



0800 31 61 81
www.alliedmedical.co.nz