

Ki Mobility Rogue2 Express Script



Customer Details					
Client					
Assessor					
Organisation					
	ACC		MOH		Other
User Weight					

Ki Mobility Rogue2 Express Script

If you only require specifications that fall under the Standard Express Package, you do not need to fill in the reverse page.

		Seat Width		
		16"	18"	20"
Seat Depth	16"	<input type="checkbox"/>	<input type="checkbox"/>	
	17"	<input type="checkbox"/>	<input type="checkbox"/>	
	18"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Express Codes

- KIROGUE21616STD
- KIROGUE21818STD
- KIROGUE21618STD
- KIROGUE22018STD
- KIROGUE21816STD

Add ROGUEEXPPACKAGE to all orders

Use KIROGUEXXXXSTD where available

Express Package: (Select optional specifications and alternatives on reverse)

- Standard Seat Sling
- Omit Transit
- 2in Depth Adjustment
- 85 degree Front Frame Bend
- Sling position 0°
- 1in Tapered Frame
- Height adjustable Medium straight back posts with no Push Handle
- Back Angle set to Standard
- Height adjustable Rigidising Bar
- Tension adjust. back upholstery
- Standard Fork
- Angle Adjustable Flip Under Footplates
- 24in Maxx Spoke Performance Pneumatic Black Wheels
- 5x1.5 Soft Roll Aluminium Castors
- Ergo Lite Handrims omit thumb guard
- Rear Wheel Camber set to 2 degrees
- COG set to 0.5"
- Caster Optimised
- Rear Wheel spacing set to 1.5" to accommodate T-Arms
- Quick Release Axles
- Push to Lock Flush Wheel Locks
- T-Arms with Standard Desk Armpads
- Rear Anti Tippers
- Iridescent Navy Frame

Please note: Express Package includes some upcharge products.

Choose Your Specifications (No additional codes are required on the below parts)

Backrest Height	Rear Seat Height																																						
<table style="width: 100%;"> <tr> <td><input type="checkbox"/> 13"</td> <td><input type="checkbox"/> 14"</td> <td rowspan="2" style="text-align: center; font-size: small;">(For 13" - 14" use .ULROGBUPTA1 if back upholstery required*)</td> </tr> <tr> <td><input type="checkbox"/> 15"</td> <td><input type="checkbox"/> 16"</td> </tr> <tr> <td colspan="3" style="text-align: center;"><input type="checkbox"/> Refer to box #1</td> </tr> </table>	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	(For 13" - 14" use .ULROGBUPTA1 if back upholstery required*)	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> Refer to box #1			<table style="width: 100%;"> <tr> <td>16W:</td> <td><input type="checkbox"/> 15"</td> <td><input type="checkbox"/> 15.5"</td> <td><input type="checkbox"/> 16"</td> <td><input type="checkbox"/> 16.5"</td> </tr> <tr> <td></td> <td><input type="checkbox"/> 17"</td> <td></td> <td></td> <td></td> </tr> <tr> <td>18W:</td> <td><input type="checkbox"/> 15"</td> <td><input type="checkbox"/> 15.5"</td> <td><input type="checkbox"/> 16"</td> <td><input type="checkbox"/> 16.5"</td> </tr> <tr> <td></td> <td><input type="checkbox"/> 17"</td> <td><input type="checkbox"/> 17.5"</td> <td><input type="checkbox"/> 18"</td> <td></td> </tr> <tr> <td>20W:</td> <td><input type="checkbox"/> 17"</td> <td><input type="checkbox"/> 17.5"</td> <td><input type="checkbox"/> 18"</td> <td><input type="checkbox"/> 18.5"</td> </tr> <tr> <td></td> <td><input type="checkbox"/> 19"</td> <td><input type="checkbox"/> 19.5"</td> <td><input type="checkbox"/> 20"</td> <td></td> </tr> </table>	16W:	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"		<input type="checkbox"/> 17"				18W:	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"		<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"		20W:	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"		<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	
<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	(For 13" - 14" use .ULROGBUPTA1 if back upholstery required*)																																					
<input type="checkbox"/> 15"	<input type="checkbox"/> 16"																																						
<input type="checkbox"/> Refer to box #1																																							
16W:	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"																																			
	<input type="checkbox"/> 17"																																						
18W:	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"																																			
	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"																																				
20W:	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"																																			
	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"																																				
<h4 style="text-align: center; background-color: #4CAF50; color: white; padding: 2px;">Backrest Angle</h4> <table style="width: 100%;"> <tr> <td><input type="checkbox"/> Std</td> <td><input type="checkbox"/> +2°</td> <td><input type="checkbox"/> +4°</td> <td><input type="checkbox"/> +6°</td> </tr> <tr> <td><input type="checkbox"/> +8°</td> <td><input type="checkbox"/> +10°</td> <td><input type="checkbox"/> -2°</td> <td><input type="checkbox"/> -4°</td> </tr> <tr> <td><input type="checkbox"/> -6°</td> <td><input type="checkbox"/> -8°</td> <td><input type="checkbox"/> -10°</td> <td></td> </tr> </table> <p>Relaxed Back <input type="checkbox"/> 0° <input type="checkbox"/> +4° <input type="checkbox"/> +9°</p>	<input type="checkbox"/> Std	<input type="checkbox"/> +2°	<input type="checkbox"/> +4°	<input type="checkbox"/> +6°	<input type="checkbox"/> +8°	<input type="checkbox"/> +10°	<input type="checkbox"/> -2°	<input type="checkbox"/> -4°	<input type="checkbox"/> -6°	<input type="checkbox"/> -8°	<input type="checkbox"/> -10°		<h4 style="text-align: center; background-color: #4CAF50; color: white; padding: 2px;">Seat to Footrest Length</h4> <table style="width: 100%;"> <tr> <td>16W:</td> <td><input type="checkbox"/> 13.5"</td> <td><input type="checkbox"/> 14"</td> <td><input type="checkbox"/> 14.5"</td> <td><input type="checkbox"/> 15"</td> </tr> <tr> <td>18W:</td> <td><input type="checkbox"/> 14"</td> <td><input type="checkbox"/> 14.5"</td> <td><input type="checkbox"/> 15"</td> <td><input type="checkbox"/> 15.5"</td> </tr> <tr> <td></td> <td><input type="checkbox"/> 16"</td> <td></td> <td></td> <td></td> </tr> <tr> <td>20W:</td> <td><input type="checkbox"/> 15"</td> <td><input type="checkbox"/> 15.5"</td> <td><input type="checkbox"/> 16"</td> <td><input type="checkbox"/> 16.5"</td> </tr> <tr> <td></td> <td><input type="checkbox"/> 17"</td> <td><input type="checkbox"/> 17.5"</td> <td></td> <td></td> </tr> </table>	16W:	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	18W:	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"		<input type="checkbox"/> 16"				20W:	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"		<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"			
<input type="checkbox"/> Std	<input type="checkbox"/> +2°	<input type="checkbox"/> +4°	<input type="checkbox"/> +6°																																				
<input type="checkbox"/> +8°	<input type="checkbox"/> +10°	<input type="checkbox"/> -2°	<input type="checkbox"/> -4°																																				
<input type="checkbox"/> -6°	<input type="checkbox"/> -8°	<input type="checkbox"/> -10°																																					
16W:	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"																																			
18W:	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"																																			
	<input type="checkbox"/> 16"																																						
20W:	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"																																			
	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"																																					
<h4 style="text-align: center; background-color: #4CAF50; color: white; padding: 2px;">Front Seat Height</h4> <table style="width: 100%;"> <tr> <td>16W:</td> <td><input type="checkbox"/> 17.5"</td> <td><input type="checkbox"/> 18"</td> <td><input type="checkbox"/> 18.5"</td> </tr> <tr> <td>18W:</td> <td><input type="checkbox"/> 18.5"</td> <td><input type="checkbox"/> 19"</td> <td><input type="checkbox"/> 19.5"</td> </tr> <tr> <td>20W:</td> <td><input type="checkbox"/> 19.5"</td> <td><input type="checkbox"/> 20"</td> <td><input type="checkbox"/> 20.5"</td> </tr> </table> <p style="text-align: center;"><input type="checkbox"/> Refer to box #7</p>	16W:	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	18W:	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	20W:	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"	<p style="text-align: center;"><input type="checkbox"/> Refer to box #4</p>																										
16W:	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"																																				
18W:	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"																																				
20W:	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"																																				

Optional Specifications

Backrest Height 1

Swap to Tall (use .ULROGBOPTL) Only use .ULROGBUPTA1 - if upholstery required*

17"
 18"
 19"
 20"

Backrest Upholstery 2

Omit Back Upholstery (add narrative)

Backrest Options 3

Std Push Handle Med (.ULROGBCKSTDMED)
 Std Push Handle Tall (.ULROGBCKSTDALL)
 Fold Down Push Handle Med (.ULROGBCKFLDMED)
 Fold Down Push Handle Tall (.ULROGBCKFLDTALL)

Footrest Options 4

High Mount Alum. Ang. Adj. Footplates A

	5"	6"	7"	8"	9"	10"	11"	12"	13"	14"
16W										
18W										
20W										

(.ULROGFPLHM1*) Based on Chair Width

Alum. Ang. Adj. Footplates B

	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"
16W									
18W									
20W									

(.ULROGFPLAL1*) Based on Chair Width

Side Guards (not available with T-Arms) 5

Clamp On Side Guards (.ULROGSDGDC)

Camber 6

0° 4° Rogue2
 (.ULROG2CAM0) (.ULROG2CAM4)

0° 4° Rogue
 (.ULROGCAM0) (.ULROGCAM4)

Castor Options 7

4 x 1.5 Soft Roll Castors (.ULROGCASSR4)
 For 16" Width, Front Seat Height 17"
 For 18" Width, Front Seat Height 18"
 For 20" Width, Front Seat Height 19"

Push Rim Options 8

Plastic Coated Push Rim (.ULEXPHRMPC)
 Std Aluminium Push Rim (.ULEXPHRMALI)
 Ergo Lite Regular (.ULEXPHRMELREG)
 Short Tabs (.ULEXPHRMSTAB)

Rear Wheel Kit - Low profile tyres 9

Rear Wheel Assembly (KIEXPWHL02)
 Includes the following:
 Spinegy spox black 24"
 Low profile - Schwalbe Marathon Plus Evolution 1"
 Ergo Lite omit thumb grip

OR

24" Maxx Spoke Low Profile Schwalbe Marathon Plus Ergo Lite Omit Thumb Grip (.ULEXPMLPFEL)
 Plastic Coated Push Rim (.ULEXPHRMPC)
 Std Aluminum Push Rim (.ULEXPHRMALI)
 Ergo Lite Regular (.ULEXPHRMELREG)
 Short Tabs (.ULEXPHRMSTAB)

Wheel Locks 10

Low Profile Scissor Locks (.ULROGWLKPLS)
 Push To Lock (.ULROGWLKPTL)

Rogue 2 Options (upcharge) 11

User Activated Anti Tips (No more than 2" dump .ULROG2UNTTP)

Accessories 12

Clear Spoke Guards (.ULROGSPGCL24)
 Impack Guards (.ULROGIMG)

Bodypoint Belt Options:

Bodypoint 2pt belt – Centre Pull Flat Mount (BPHB205-M46-A1) + (BPHW200-90-2)
 Bodypoint 4pt belt Centre Pull Flat Mount (BPHB405-M46-A1) + (BPHW200-90-2)
 Bodypoint Evoflex Pelvic Stabiliser (BPEB205-M46) + (BPHW200-90-2)
 Bodypoint Aeromesh Calf Strap (BPSP103M)

Notes

