



**ALLIED
MEDICAL
LIMITED**

Knee Bracing & Support Solutions



	UNIVERSAL PATELLA SUPPORT	AML SUPER KNEE SPLINT	X ACTION KNEE BRACE	ACTIVE X PATELLA	AML KNEE R.O.M.	WRAPAROUND HINGED KNEE SUPPORT	LATERAL PATELLA STABILISER	KNEE RANGER	REHAB AIRLIGHT	CHECK BRACE	UNLOADER ONE OTS KNEE BRACE
BRIEF DESCRIPTION	The patella support is designed to hold the patella in position. Worn beneath patella, it is suitable for the relief of Osgood-Schlatter syndrome, chondromalacia and tendonitis.	A comfortable and sturdy knee immobiliser for post surgery or after knee injury.	This is an economical elastic knee support with a crossover back. Provides extra support around the knee.	The Allard Active X Patella is designed to stabilise the patella and prevent dislocation and subluxation. With a supportive built-in pad, the knee is kept in place while moving.	The AML Knee R.O.M. allows for post operative / injury immobilisation and protected range of motion.	A low profile brace with Dual Action Hinges moulded into the brace and adjustable strapping to give extra support.	This pull on sleeve has a sewn in 'J' buttress with dynamic side pull compression strap to support patella malalignment, dislocation and subluxation.	This brace offers post-operative immobilisation with protected range of motion.	The Rehab Airlight is ideal as a post operative or post injury brace when ROM control and immobilisation is required.	Comfortable knee brace that controls hyper-extension of the knee. The brace helps with stability and support when moving and when resting.	The flexible, perforated shells conform easily to individual anatomy, giving the Unloader One OTS an excellent out-of-the-box fit.
CONSTRUCTION	Durable vinyl/tricot laminate with velcro closure	Foam and mesh with velcro	Elastic/cotton	Ventilated Active X-material	Foam leg wraps with adjustable velcro straps	Neoprene lined with lycra	Neoprene	Foam	Air light padding, foam	Synthetic leather laminated with polyurethane foam and brushed polyester tricot	Ossur Sensil Silicone technology and aerospace
COLOUR	Beige	White/Blue	Beige	Black	Black	Black	Black	Black	Black	Black	Grey
SIZES	Universal	Small - X Large	Small - X Large	X Small - Large	Short or Long	Small - XXX Large	Small - X Large	Short or Long	Universal	Universal	Small - Large
ORIENTATION	Unilateral	Unilateral	Unilateral	Unilateral	Unilateral	Unilateral	Left or Right	Unilateral	Left or Right	Unilateral	Left or Right
LENGTH	4.3cm	30cm	20cm	14cm	55cm or 66cm	30cm	32.5cm	42.5cm or 65cm	49 - 60cm	45cm	44cm
PRODUCT CODE	PR80300	AML1320x	AMLXAKx	AL28055001x	AMLROMKx	OS21200x	PR944xx	PR942x0	OSB220x00000	AL286520000	OSB240xxxxxx
BRAND	DJ Orthotics Procure	DJ Orthotics	Allied Medical	Allard International	Allied Medical	Ossur	DJ Orthotics Procure	DJ Orthotics Procure	Ossur	Allard International	Ossur

Information correct as of February 2019. From time to time, product specifications are changed by manufacturers without prior notice

Allied Medical Limited PH: 0800 31 61 81 sales@alliedmedical.co.nz