



**ALLIED  
MEDICAL  
LIMITED**

## Leckey MyWay+ and MyWay Pedal

<b>Name:</b>	
<b>Date:</b>	
<b>Email:</b>	
<b>Contact Number:</b>	
<b>Invoice Address:</b>	
<b>Delivery Address:</b>	

**Purchase Order**

Quote

Order



Item Code	Description	Quantity	Available for Trial
<b>Frame</b>			
<input type="checkbox"/>	LK209-1600 MyWay+ S1 MyWay Frame		✓
<input type="checkbox"/>	LK209-2600 MyWay+ S2 MyWay Frame		✓
<input type="checkbox"/>	LK209-3600 MyWay+ S3 MyWay Frame		✓
<b>Harness</b>			
<input type="checkbox"/>	LK173-4100-04 MyWay S1 Purple MyWay Harness		✓
<input type="checkbox"/>	LK173-4200-07 MyWay S2 Green MyWay Harness		✓
<input type="checkbox"/>	LK173-4300-02 MyWay S3 Blue MyWay Harness		✓
<input type="checkbox"/>	LK173-4400-09 MyWay S4 Red MyWay Harness		✓
<input type="checkbox"/>	LK173-4500-06 MyWay S5 Orange MyWay Harness		✓
<input type="checkbox"/>	LK173-4600-08 MyWay S6 Grey MyWay Harness		✓
<input type="checkbox"/>	LK173-1800 MyWay Hoist Straps (one size)		✓
<b>Headrest</b>			
<input type="checkbox"/>	LK173-1680 MyWay Headrest		✓
<input type="checkbox"/>	LK209-1680 MyWay+ Complex Headrest		
<b>User Handles</b>			
<input type="checkbox"/>	LK209-1667 MyWay+ Size 1 User Handles (Pair)		✓
<input type="checkbox"/>	LK209-2667 MyWay+ Size 2 User Handles (Pair)		✓
<input type="checkbox"/>	LK209-3667 MyWay+ Size 3 User Handles (Pair)		✓
<input type="checkbox"/>	LK209-1740 MyWay+ Size 1 Complex User Handles incl. Ball Grip (Pair)		
<input type="checkbox"/>	LK209-2740 MyWay+ Size 2 Complex User Handles incl. Ball Grip (Pair)		
<input type="checkbox"/>	LK209-3740 MyWay+ Size 3 Complex User Handles incl. Ball Grip (Pair)		
<input type="checkbox"/>	LK209-1319 MyWay+ Cylinder Grip (Pair)		
<b>Adult Handles</b>			
<input type="checkbox"/>	LK209-1691 MyWay+ S1 Adult Handles		✓
<input type="checkbox"/>	LK209-2691 MyWay+ S2 Adult Handles		✓
<b>Ankle Guides</b>			
<input type="checkbox"/>	LK209-1700 MyWay+ S1 Ankle Guides (L&R)		✓
<input type="checkbox"/>	LK209-2700 MyWay+ S2 Ankle Guides (L&R)		✓
<input type="checkbox"/>	LK209-3700 MyWay+ S3 Ankle Guides (L&R)		✓

Saddle				
				
<input type="checkbox"/>	LK173-1671	MyWay S1 Saddle (recommended for MyWay S1 & S2 frames)		✓
<input type="checkbox"/>	LK173-2671	MyWay S2 Saddle (recommended for MyWay S2 & S3 frames)		✓
<input type="checkbox"/>	LK209-1671	MyWay+ S1 Narrow Saddle (recommended for MyWay S1 & S2 frames)		
<input type="checkbox"/>	LK209-2671	MyWay+ S2 Narrow Saddle (recommended for MyWay S2 & S3 frames)		
<input type="checkbox"/>	LK209-1679	MyWay+ S1 Drop (fits with MyWay S1 saddle only)		
<input type="checkbox"/>	LK209-2679	MyWay+ S2 Drop (fits with MyWay S2 saddle only)		
<input type="checkbox"/>	LK209-3679	MyWay+ S1 Drop (fits with MyWay S2 saddle only)		
MyWay Pedal (Size 2 and 3 MyWay only)				
				
<input type="checkbox"/>	LK173-1750	Smart Pedal Walker (one size) requires sandals		
<input type="checkbox"/>	LK152-1600	Sandals Size 1		✓
<input type="checkbox"/>	LK152-2600	Sandals Size 2		✓
<input type="checkbox"/>	LK152-3600	Sandals Size 3		✓
<input type="checkbox"/>	LK173-2305	Small Ankle Wraps (use with sandals)		
<input type="checkbox"/>	LK173-3305	Large Ankle Wraps (use with sandals)		

**Important Note - MyWay accessories not shown on this order form are no longer compatible or available to order with MyWay+**

Frame Size	1	2	3			
Age [approx.]	1-5	4-11	9-16			
User Shoulder Height [to floor] (mm)	590-870	870-1190	1070-1490			
User Saddle Height [inside leg] (mm)	240-410	400-610	550-810			
Max. User Weight (kg)	25	50	80			
Frame Weight (kg)	14	7.5	11			
Upper Support Unit Weight (kg)	Size 1 frame & upper support unit do not separate	10	11			
Frame Width (mm)	570	670	720			
Frame Length (mm)	840	930	1180			
Harness Size	1	2	3	4	5	6
Age [approx.]	1-4	2-6	5-9	8-12	12-15	Up to 16
Waist at Navel Level (cm)	46-56	50-60	55-65	60-70	70-80	80-90
Max. User Weight (kg)	50	50	50	80	80	80
Colour	Purple	Green	Blue	Red	Orange	Grey

Comments