



LECKEY

MyWay+

Adventure in every step

Supporting Children

Introducing the MyWay+

MyWay+ is an open frame walker which cleverly supports children in an upright position, improving stepping and maximising opportunities for exploration, interaction and participation.

A bespoke frame combined with a wraparound hoistable harness and a wide range of accessories, ensures the MyWay+ can be tailored to meet every need and ability.

Scan with your camera to view the new MyWay+ in action



Support To Go Your Own Way

MyWay+ continues to lead the way when it comes to opening a world of independence and exploration for children.

Building on the very best of the MyWay walking frame, we are excited to introduce MyWay+, reimagined with a wide range of new features and accessories to suit every need. By allowing the MyWay+ to be tailored to each child's individual needs, we have opened up even more possibilities for children to play, have fun and explore without limitations.

Early intervention is key to optimising a child's potential and therefore we have resized our size 1 frame to better accommodate younger children's body proportions and really give them the best opportunity to develop their physical, cognitive, emotional and social skills.



“
Emilia absolutely loves the freedom that her MyWay+ walker offers her and being able to join in with her siblings when they are playing together. The MyWay+ keeps all her joints and muscles mobile and gives her the best opportunities for movement and independence.

Emilia's mum



Tailored Approach To Every Child

Available in 3 frame sizes with 6 colour-coded harnesses, MyWay+ brings individuality to the world of walking frames. Each frame comes with our robust new indoor/outdoor wheels ensuring there's no limit to which path your little one can take. These sturdy wheels with *side-step technology come with anti-roll back, direction lock and variable resistance as standard creating that individual set-up to get everyone moving.

A wide range of accessories can be combined to create a bespoke frame, tailored to every child's needs and abilities.



Learn about the MyWay+ accessories



OPTIONS INCLUDE:

- 1 Castors – Indoor and Outdoor (Robust, anti-roll back, direction lock (including *sideways lock), varied resistance) – to suit a child's individual needs



Learn about the MyWay+ Castors

- 2 Headrests (complex and standard) – to act as a prompt or offer more intimate occipital support.
- 3 User Handles (complex and standard) – multi-adjustable support for the upper limbs may aid stepping
- 4 Saddle drop – proximal support to help reduce hip adduction
- 5 Saddle (standard and slim) – to aid with transfers of provide additional off-loading of lower limbs

- 6 Three options of frame size
- 7 Tailored fit from toddler to teen for frames and harnesses
- 8 Overall and fine height adjustment of the spine
- 9 Detachable and height adjustable adult handles – to enable parents or carers to safely push and steer the frame, particularly useful with the pedal accessory
- 10 Ankle guides – to guide lower limb movement away from scissoring or asymmetric stepping pattern

*Patent Pending - United Kingdom Application No 2302760.0



Unique Corset Like Harness

The revolutionary corset-like harness on MyWay+ has been designed to provide 360° support to promote an optimal stepping pattern as well as being hoistable. By focusing on function and a more natural upright movement pattern, the MyWay+ provides all the support and reassurance children need to do things their way - from kicking a football or participating in playground games, to walking hand-in-hand with their siblings.

- Hug-like effect increases sensory feedback, dampens involuntary movement and provides comfort, support and reassurance.
- Provides added core stability which in turn can improve lower limb movements and head control
- Hoistable harness encourages easier transitions in and out of the frame



Meet Harry



Open Frame Design

Promoting freedom and independence are at the heart of MyWay+. That's why we have a minimalist, open-frame design that provides support without restriction and opens up a world of possibilities.

- Open, easy-access frame encourages interaction and participation, whether at home, in the playground or in the classroom.
- Newly developed indoor/outdoor castors enable children to be active in lots of different environments, encouraging children of all abilities to meet the recommended physical activity guidelines.



Meet Tom





Supporting Therapy Practice

Children with additional needs often have weaker muscles, especially around their hips, and this can impact on their gait, function and also hip joint development. The new *side-step function of the MyWay+, locks the wheels of the frame in a sideways position providing opportunity to target and strengthen hip abductor muscles whilst the child is fully supported in the frame.



Scan with your camera to see Willow's fantastic improvement in hip muscle strength after using the *side-stepping function

Using the MyWay Pedal accessory enables children to target and strengthen their anti-gravity muscles (calf, quads and gluteals) in an upright, weight-bearing position - unlike a trike which encourages a flexed position and overuse of the hip flexors and movement in a circular pattern. The MyWay Pedal provides a sensori-motor experience more similar to walking as it facilitates an elliptical movement pattern, whilst encouraging improved lower limb alignment and symmetry. It is suitable for children with GMFCS III to V and can be used as part of a child's walking therapy programme or trialed as a precursor to stepping.



Learn about the MyWay Pedal

There's always more to MyWay+

Multi-functional castors, lightweight frame, prone, gross and fine height adjustment combine to help initiate stepping and control movement.

Frame Size	1	2	3
Age (approx.)	1-5	4-11	9-16
User Shoulder Height [to floor] (mm)	590-870	870-1190	1070-1490
User Saddle Height [inside leg] (mm)	240-410	400-610	550-810
Max. User Weight (kg)	25	50	80
Frame Weight (kg)	14	7.5	11
Upper Support Unit Weight (kg)	Size 1 frame & upper support unit do not separate	10	11
Frame Width (mm)	570	670	720
Frame Length (mm)	840	930	1180

Harness Size	1	2	3	4	5	6
Age (approx.)	1-4	2-6	5-9	8-12	12-15	Up to 16
Waist at Navel Level (cm)	46-56	50-60	55-65	60-70	70-80	80-90
Max. User Weight (kg)	50	50	50	80	80	80
Colour	Purple	Green	Blue	Red	Orange	Grey

To arrange a product demonstration/assessment or for more information, please contact us:



LECKEY