



MyWay+ Measurement Guide


Child name: _____

Date: _____

Location: _____

Therapist: _____

Notes: _____



- A** Inside leg – heel: _____
- B** Top of shoulder – heel: _____
- C** Waist circumference at naval level (for harness size):

Frame Size	1	2	3
Age (approx.)	1-5	4-11	9-16
User Shoulder Height (to floor) (mm)	590-870	870-1190	1070-1490
User Saddle Height (inside leg) (mm)	240-410	400-610	550-810
Max. User Weight (kg)	25	50	80
Frame Weight (kg)	14	7.5	11
Upper Support Unit Weight (kg)	Size 1 frame & upper support unit do not separate		11
Frame Width (mm)	570	670	720
Frame Length (mm)	840	930	1180

Harness Size	1	2	3	4	5	6
Age (approx.)	1-4	2-6	5-9	8-12	12-15	Up to 16
Waist at Navel Level (cm)	46-56	50-60	55-65	60-70	70-80	80-90
Max. User Weight (kg)	50	50	50	80	80	80
Colour	Purple	Green	Blue	Red	Orange	Grey