






Leckey Squiggles Seat

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order

Purchase Order	
----------------	--



	Item Code	Description	Quantity	Available for Trial
Seat				
				
<input type="checkbox"/>	120-610-MIN	Squiggles Seat Squiggles Shell-Configured in Minimum Seat Depth Set-up		
Interface				
				
<input type="checkbox"/>	194-600	Squiggles Seat Interface for Sunrise bases		
High Low Base				
				
<input type="checkbox"/>	120-800	Squiggles Seat Hi low chassis - foot pedal		
<input type="checkbox"/>	122-710	Mygo/Squiggles Standard Black Push handle		✓
<input type="checkbox"/>	122-2710	Mygo/Squiggles Height Adjustable Push Handle (Grey Metal)		

Mygo Base



Indoor



Indoor/Outdoor

<input type="checkbox"/>	120-818	Interface to Attach Squiggles Seat to Mygo Chassis		
<input type="checkbox"/>	102-2753-G	Hi-Low Mygo Grey Indoor Manual Chassis - With Height Adjustable Handle		
<input type="checkbox"/>	102-2754-G	Hi-Low Mygo Grey Indoor Powered Chassis - With Height Adjustable Handle		
<input type="checkbox"/>	102-2766-G	Hi-Low Mygo Grey Indoor Hydraulic Chassis - With Height Adjustable Handle		
<input type="checkbox"/>	102-2767-G	Hi-Low Mygo Grey Indoor/Outdoor Manual Chassis - With Height Adjustable Handle		
<input type="checkbox"/>	102-2768-G	Hi-Low Mygo Grey Indoor/Outdoor Hydraulic Chassis - With Height Adjustable Handle		
<input type="checkbox"/>	102-2769-G	Hi-Low Mygo Grey Indoor/Outdoor Powered Chassis - With Height Adjustable Handle		
<input type="checkbox"/>	102-2753	Hi-Low Mygo Black Indoor Manual Chassis - With Push Handle		
<input type="checkbox"/>	102-2754	Hi-Low Mygo Black Indoor Powered Chassis - With Push Handle		
<input type="checkbox"/>	102-2766	Hi-Low Mygo Black Indoor Hydraulic Chassis - With Push Handle		
<input type="checkbox"/>	102-2767	Hi-Low Mygo Black Indoor/Outdoor Manual Chassis - With Push Handle		
<input type="checkbox"/>	102-2768	Hi-Low Mygo Black Indoor/Outdoor Hydraulic Chassis - With Push Handle		
<input type="checkbox"/>	102-2769	Hi-Low Mygo Black Indoor/Outdoor Powered Chassis - With Push Handle		

Tray



Tray




Padded Tray Cover



Grab Rail

<input type="checkbox"/>	120-672	Squiggles Seat tray		
<input type="checkbox"/>	120-688	Squiggles Seat tray 50 mm extended height		
<input type="checkbox"/>	117-769	MYGO/Squiggles Grab Rails		
<input type="checkbox"/>	188-1854-05	Squiggles Black Padded Tray Cover		

Covers



Blue

Green

Pink

<input type="checkbox"/>	120-710-04	Squiggles Seat covers - pink (Includes femoral seat and backrest)		
<input type="checkbox"/>	120-710-03	Squiggles Seat covers - blue (Includes femoral seat and backrest)		
<input type="checkbox"/>	120-710-01	Squiggles Seat covers - green (Includes femoral seat and backrest)		

Shoulder Support



<input type="checkbox"/>	120-702	Squiggles Seat Shoulder Support Assembly		
<input type="checkbox"/>	120-828-04	Squiggles Seat Pink Shoulder Cushion		
<input type="checkbox"/>	120-828-03	Squiggles Seat Blue Shoulder Cushion		
<input type="checkbox"/>	120-828-01	Squiggles Seat Green Shoulder Cushion		

Leg Guide				
				
<input type="checkbox"/>	120-856-04	Squiggles Seat Leg Guide Pad Assembly - pink		
<input type="checkbox"/>	120-856-03	Squiggles Seat Leg Guide Pad Assembly - blue		
<input type="checkbox"/>	120-856-01	Squiggles Seat Leg Guide Pad Assembly - green		✓
Footrest and Pommel				
 				
<input type="checkbox"/>	120-641	Squiggles Seat pommel		
<input type="checkbox"/>	120-639	Squiggles Seat Footrest		
Back Spacer				
				
<input type="checkbox"/>	120-779-04	Squiggles Seat 1" back spacer pad & flipaway lateral support - Pink		
<input type="checkbox"/>	120-779-03	Squiggles Seat 1" back spacer pad & flipaway lateral support - Blue		
<input type="checkbox"/>	120-779-01	Squiggles Seat 1" back spacer pad & flipaway lateral support - Green		
<input type="checkbox"/>	120-777-04	Squiggles Seat 1" back spacer pad & rigid lateral support - Pink		
<input type="checkbox"/>	120-777-03	Squiggles Seat 1" back spacer pad & rigid lateral support - Blue		
<input type="checkbox"/>	120-777-01	Squiggles Seat 1" back spacer pad & rigid lateral support - Green		
Chest Support				
 				
<input type="checkbox"/>	137-808-05	Oval Chest Harness - Black		
<input type="checkbox"/>	137-808-05VIN	Oval Chest Harness - Black Vinyl		
<input type="checkbox"/>	155-N001-L600	Small Butterfly Harness & attachment straps		✓
<input type="checkbox"/>	155-N001-L601	Medium Butterfly Harness & attachment straps		✓
Foot Support				
   				
	Sandals	Booties	Ankle Huggers	Ankle Wraps
<input type="checkbox"/>	152-1600	Sandals Size 1		✓
<input type="checkbox"/>	150-1852	Extra Small Booties		
<input type="checkbox"/>	AH01	Ankle Huggers		
<input type="checkbox"/>	173-2305	Small Ankle Wraps (use with sandals)		

<input type="checkbox"/>	120-816	1Inch Footplate Raiser for Single Footplate Head Support		✓
				
	Contoured	Contoured Occipital	Flat	Laterals (for use with flat)
<input type="checkbox"/>	137-699-05	Contoured headrest Assembly and Cushion		
<input type="checkbox"/>	137-699-05VIN	Contoured headrest Assembly and Cushion - Vinyl		
<input type="checkbox"/>	117-699-05	Contoured Occipital Headrest		
<input type="checkbox"/>	117-699-05VIN	Contoured Occipital Headrest - Vinyl		
<input type="checkbox"/>	117-851-05	Flat Headrest Assembly and Cushion		✓
<input type="checkbox"/>	117-851-05VIN	Flat Headrest Assembly and Cushion - Vinyl		
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with covers		
<input type="checkbox"/>	120-893-05VIN	Flat Headrest Laterals with covers Vinyl		
Pelvic Support				
				
	Pelvic Harness	Pelvic Spacer		
<input type="checkbox"/>	120-844-03	Squiggles pelvic harness small - blue		
<input type="checkbox"/>	120-844-01	Squiggles pelvic harness small - green		
<input type="checkbox"/>	120-844-04	Squiggles pelvic harness small - pink		
<input type="checkbox"/>	117-848-03	Squiggles pelvic harness medium - blue		
<input type="checkbox"/>	117-848-01	Squiggles pelvic harness medium - green		✓
<input type="checkbox"/>	117-848-04	Squiggles pelvic harness medium - pink		
<input type="checkbox"/>	120-789	MYGO/Squiggles Pelvic Harness spacer pads (pair)		
Laterals				
				
	Rigid Laterals	Flipaway Laterals	Flipaway Laterals with Protraction Pads	
<input type="checkbox"/>	120-656	Squiggles Seat rigid lateral support - Hardware		
<input type="checkbox"/>	120-855	Squiggles Seat flipaway lateral support Hardware		
<input type="checkbox"/>	120-838-03	Squiggles Seat lateral support covers - blue		
<input type="checkbox"/>	120-838-01	Squiggles Seat lateral support covers - green		✓
<input type="checkbox"/>	120-838-04	Squiggles Seat lateral support covers - pink		
<input type="checkbox"/>	117-644-05	Flipaway Laterals with Protraction Pads		

Product Sizing			
		mm	Inches
Age (Approx)	1-5 Years		
Stature	Min	750	29.5
	Max	1100	43.3
Max User Weight (kg/Lbs)		22Kg	48.4Lbs
Seat Depth (Backrest to back of knee)	Min	190	7.5
	Max	290	11.4
Seat width	Min	160	6.3
	Max	260	10.2
Hip Width (Distance between hip guides)	Min	160	6.3
	Max	220	8.7
Backrest Height (Sitting shoulder height)	Min	320	12.6
	Max	380	15
Backrest Height (with shoulder support)	Min	420	16.5
	Max	480	18.9
Chest Width (Distance between laterals)	Min	160	6.3
	Max	90	3.5
Distance between pommel & side of seat	Min	190	7.5
Seat to Footplate (Sole of foot to back of knee)	Min	122.5	4.8
	Max	650	25
Seat to Floor (Top of Seat to Floor (Hi-Low Chassis))	Min	270	10
	Max	650	25
Backrest Angle	Prone	10°	10°
	Recline	25°	25°
Footplate Angle	Plantarflexion	10°	10°
	Dorsiflexion	10°	10°
Tray Size		460/560	18.1/22

Comments