



**Leckey  
Totstander**

Order no.:

-----

Provide estimate only     Order

Attachment for custom fabrication enclosed

Billing address/customer number

Delivery address

Sales person

Name

.....  
Date

Signature

# Leckey Totstander



## Standard Equipment

- Stable frame of steel and wood elements
- Adjustable chest and knee supports, lap belt
- Adjustable sandals

Standard model		
LKLFS/T-MW	<input type="checkbox"/>	<b>Stander for children from approx. 1 to 5 years</b> , cover colour: Meadow for vertical standing

Technical data		
<b>Age</b> (years)	<b>1 - 5</b>	
<b>Trunk belt height</b>	<b>Min.</b>	<b>470/18.5</b>
	<b>Max.</b>	<b>640/26</b>
<b>Base plate</b>	<b>Depth</b>	<b>450/18</b>
	<b>Width</b>	<b>450/18</b>