



Company: _____

Therapist: _____ Date: _____

Acc MoH

Marked For: _____

Liberty FT E1161 _____ \$
113 kg weight capacity. Tilt range: 0° - 20°.

Ship To: _____

Address: _____

City: _____ State: _____

Phone: _____

Email/Fax: _____

Notes: _____ Effective 21/04/2021

Frame

Seat Width Range / Seat Width

Select Seat Width within desired Seat Width Range

Seat Width Range	Seat Width				
Small 14" - 18"	<input type="checkbox"/> 14 N/C	<input type="checkbox"/> 15 N/C	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C
Medium 16" - 20"	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 \$ E2201
Large 18" - 22"	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 \$ E2201	<input type="checkbox"/> 21 \$ E2201	<input type="checkbox"/> 22 \$ E2201

Assistive Gas Strut Set-Up

Client Weight	Setting	Spring Type (Newtons)
<input type="checkbox"/> > 90 kgs	1	450 N
<input type="checkbox"/> 79 - 90 kgs	2	450 N
<input type="checkbox"/> 68 - 79 kgs ¹	3	450 N
<input type="checkbox"/> 68 - 79 kgs ²	1	300 N
<input type="checkbox"/> 56 - 68 kgs	2	300 N
<input type="checkbox"/> < 56 kgs	3	300 N

¹ Standard on chairs less than 16" deep

² Standard on chairs 16" deep or larger

Seat Depth Range / Seat Depth

Select Seat Depth within desired Seat Depth Range

Seat Depth Range	Seat Depth		
Short 14" - 16"	<input type="checkbox"/> 14 N/C	<input type="checkbox"/> 15 N/C	<input type="checkbox"/> 16 N/C
Medium 16" - 18"	<input type="checkbox"/> 16 N/C	<input type="checkbox"/> 17 N/C	<input type="checkbox"/> 18 N/C
Long 18" - 20"	<input type="checkbox"/> 18 N/C	<input type="checkbox"/> 19 N/C	<input type="checkbox"/> 20 \$ E2203

Seat Option

- Standard Seat Sling _____ Std
- Omit Seat _____ N/C

Transit Option (Tie Downs)

- Transit System K0108 _____ \$

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post and Height Range

Back Post	Price	Height Range	Back Height				
Stroller Handle Back Post ¹	Std	Short	<input type="checkbox"/> 20				
		Tall	<input type="checkbox"/> 24				
Height Adjustable Straight with Push Handle	\$	Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
		Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20
		Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22
Height Adjustable 8° Bend with Push Handle ²	\$	Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
		Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20
		Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22

¹ Not available with Back Upholstery

² Not available with IV Pole

Stroller Handle Back Post Option

Available only with Stroller Handle Back Post

Removable Stroller Handle Extension _____ \$

Backrest Options

Fixed (Non-Folding) _____ \$

Angle Adjustable (Folding) _____ \$

Back Upholstery

Standard _____ Std

Tension Adjustable E2611 _____ \$
Not available on Seat Widths less than 16" combined with 21" and 22" Backrest Heights

Omit Back Upholstery _____ N/C

Footrest

Hanger

60° _____ Std
Available only with Extension Tube Hangers

70° _____ N/C

80° _____ N/C

90° _____ N/C
Available only with Front Mount Hangers. Not available with 8" Casters combined with 3" Extension

Elevating Leg Rests (ELR) E0990, Rental K0195 _____ \$

Pro Elevating Leg Rests (PRO ELR) E0990, Rental K0195 _____ \$

Omit Hangers _____ N/C

PRO ELR Cover

Gel Pad Cover _____ \$

Footplate

Composite _____ Std
Not available with Front Mount Hanger or PRO ELR

Composite Angle Adjustable _____ \$
Not available with Front Mount Hanger or 14" Seat Width

Aluminum Angle Adjustable K0040 _____ \$

Aluminum Locking Multi-Angle Adjustable K0108 _____ \$
Available only with Front Mount Hanger

One Piece Flip Up Angle Adjustable K0040 _____ \$
Not available with Front Mount Hanger, ELR or PRO ELR

Footrest continued

Footrest / Extension / Seat to Footrest Length

Select Seat to Footrest Length in row corresponding to desired Footrest and Extension. Footrest Length of at least 2.5" shorter than Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length												
Extension Tube Hangers (60°, 70°, 80°)	Short	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5					
	Medium	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5				
	Long	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18				
Front Mount Hangers ¹	Standard Selections										Requires 3" Extension			
	Short (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5						<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	
		<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8								<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11
	Medium (70°, 80°, 90°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5					<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5
		<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10						<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
	Long (70°, 80°)	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5		
<input type="checkbox"/> 9		<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15				
Pro ELR Low Pivot Setting	Short	<input type="checkbox"/> 10.5 - 15.5												
	Long	<input type="checkbox"/> 13 - 19												
Pro ELR High Pivot Setting	Short	<input type="checkbox"/> 8.5 - 17												
	Long	<input type="checkbox"/> 11.5 - 17.5												

¹ See Minimum Seat Height chart below

Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

Heel Loop

Not available with One Piece Flip Up Angle Adjustable Footplate

- Heel Loops E0951 _____ \$
- Omit Heel Loops _____ N/C

Calf Strap

- Velcro Adjustable K0038 _____ \$
- Padded Velcro Adjustable K0038 _____ \$
- Bodypoint Padded AeroMesh K0038 _____ \$

3" Extension

Available only with Front Mount Hangers. Provided as parts.

- 3" Extension for Front Mount Hanger _____ N/C

Residual Limb Support

Residual Limb Support E1020

- Left _____ \$
- Right _____ \$
- Omit Left Hanger _ N/C
- Omit Right Hanger _ N/C

Cover

- Smooth Air Mesh
- Smooth Air Mesh


Pad Size

- 5 x 8 9 x 14
- 5 x 8 9 x 14
- 7 x 10 11 x 10
- 7 x 10 11 x 10
- 9 x 10 11 x 14
- 9 x 10 11 x 14

Footrest continued

Hanger Release

- 4-Way - In line _____ Std
See reference A to the right. Not available with ELR or PRO ELR.
- 4-Way - Out _____ N/C
See reference B to the right



Seat Height, Casters and Rear Wheels

Seat Height / Caster Size / Rear Wheel Size

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size			
	4"	5"	6"	8"
12"	Seat Height			
	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16.5 ⁴
	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 17
	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 18
	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 19 ²
	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	
		<input type="checkbox"/> 18	<input type="checkbox"/> 19 ²	
16"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16.5 ⁴
	<input type="checkbox"/> 14 ¹	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 17
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 18
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 19
		<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19.5
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 20
		<input type="checkbox"/> 18	<input type="checkbox"/> 19	
		<input type="checkbox"/> 18.5		
20"	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16.5 ⁴
	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 17
	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 18
	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 19
	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19.5
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 20
		<input type="checkbox"/> 18	<input type="checkbox"/> 19	
		<input type="checkbox"/> 18.5		

Rear Wheel Size	Caster Size			
	4"	5"	6"	8"
22" ³	Seat Height			
	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16.5 ⁴
	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 17
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 18
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18.5
	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 19
		<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19.5
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 20
		<input type="checkbox"/> 18	<input type="checkbox"/> 19	
24" ³	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 16.5 ⁴
	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 17
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 18
		<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 19
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 19.5
		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 20
		<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	
		<input type="checkbox"/> 19		

¹ Available only with Low Profile tires

² Not available with Full Poly tire

³ Not available with Low Profile tires in combination with Mag wheel

⁴ Available only with 1.5" Poly casters

Seat Height, Casters and Rear Wheels continued

Caster Type

- 1" Poly _____ Std
- 1" Poly Aluminum _____ \$
Available only with 4" and 5" Caster Sizes
- 1" Pneumatic E2214, E2215 _____ \$
Available only with 6" and 8" Caster Sizes
- 1.5" Poly _____ \$
- 1.5" Soft Roll Aluminum _____ \$
Not available with 8" Caster Size
- 2" Poly _____ \$
Available only with 6" Caster Size
- 2" Pneumatic w/ Foam Insert E2217 _____ \$
Available only with 8" Caster Size

Rear Wheel Type

- 24" Spoke _____ N/C
- Spoke (20" or 22") _____ \$
- Mag _____ Std

Rear Tire Type

Full Profile

- Street Pneumatic E2211, E2212 _____ Std
- Full Poly E2218 _____ N/C
- Pneumatic with Airless Insert E2211, E2213 _____ \$
- Euro Trek Knobby E2211, E2212 _____ \$
Available only on 24" wheel size

Low Profile

- Low Poly E2218 _____ N/C
Available only on 16" wheel size
- Sentinel Solid Tire E2218 _____ \$
Not available on 12" or 16" wheel sizes or with 20" and 22" Mag wheels
- Schwalbe Marathon Plus Evolution E2211, E2212 _____ \$
Available only on 24" spoke wheel

Wheel Locks

Wheel Lock

- Push to Lock _____ Std
- Pull to Lock _____ N/C
Not available with 12" or 16" Full Poly tires
- Hemi Wheel Lock _____ \$
Extension Handle included
 - Left Right
- Drum Brake _____ \$

Armrests

Armrest Please select Arm Pad option

- Height Adjustable T-Arm (8.5" - 13.5") E0973 _____ \$
- Angle Adjustable Locking Extendable Flip Up E0973 _____ \$
Available only with 12" and 16" wheel sizes and Stroller Handle Back Post
- Omit Armrests _____ N/C
- Tubular Flip Up _____ Std
Available only with Stroller Handle Back Post. Not available with Arm Pads.

Handrim

Handrims not available on 12" or 16" wheel sizes

Standard Handrims

- Aluminum Anodized _____ Std
- Plastic Coated _____ \$
- Projection (8 Vertical Projections) _____ \$
Available only with full profile tires in combination with 22"-24" Spoke or 24" Mag wheel
- Omit Handrims _____ N/C

Natural Fit Handrims

Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.

- Natural Fit LT Regular Grip _____ \$
- Natural Fit Standard Regular Grip _____ \$
- Natural Fit LT Super Grip _____ \$
- Natural Fit Standard Super Grip _____ \$
- Natural Fit LT Omit Grip _____ \$
- Natural Fit Standard Omit Grip _____ \$

Handrim Option

Not available on Plastic Coated or Projection Handrims

- 3M Micro Grip _____ \$

Axle

- Fixed _____ N/C
- Quick Release _____ Std

- Dual Drum Brake _____ \$

Available only with 12" Mag, 16" Mag, or MAXX Performance Spoke wheels. Axle is quick release.

- Push to Lock Pull to Lock

- Omit Wheel Locks _____ N/C

Wheel Lock Options

- Extension Handles E0961 _____ \$
Available only with Push to Lock or Pull to Lock Wheel Locks
- Grade Aids E0974 _____ \$
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires

Arm Pad

- Standard Arm Pad _____ N/C Desk Full
Not available with Angle Adjustable Locking Extendable Flip Up Armrests in combination with 20"-24" wheel sizes
- Waterfall Arm Pad _____ \$ Desk Full
Not available with Angle Adjustable Locking Extendable Flip Up Armrests
- Foam Grip _____ N/C Desk Full
Available only with Angle Adjustable Locking Extendable Flip Up Armrests

Accessories

Removable Side Guard		O₂ Holder	
<input type="checkbox"/> Composite Side Guards K0108	\$	<input type="checkbox"/> O ₂ Holder E2208	\$
<i>Not available with T-Arms</i>		<i>Mounted to the left side. Not available with Fold Down Back.</i>	
Anti-Tip		Other Accessories	
<input type="checkbox"/> Rear Anti-Tips E0971	\$	<input type="checkbox"/> Tool Kit	\$
Positioning Belt		<input type="checkbox"/> Removable Underseat Pouch	\$
<input type="checkbox"/> 1 1/2" Auto Buckle E0978	\$	<input type="checkbox"/> Ki Mobility Backpack	\$
<input type="checkbox"/> 1 1/2" Padded Auto Buckle E0978	\$	<input type="checkbox"/> Black Neoprene Impact Guards	\$
<input type="checkbox"/> 2" Airline Buckle E0978	\$	<input type="checkbox"/> Clear Spoke Guards K0065	\$
Bodypoint Evoflex - Pelvic Stabilizer E0978	\$	<i>Not available with Mag wheel</i>	
Size: <input type="checkbox"/> Small (16-24") <input type="checkbox"/> Medium (19-29") <input type="checkbox"/> Large (25-37")		<input type="checkbox"/> Black Spoke Guards K0065	\$
Belt Mounting Kit		<i>Not available with 20" wheel size or Mag wheel</i>	
<input type="checkbox"/> Seat Belt Mounting Hardware	\$	<input type="checkbox"/> Cane and Crutch Holder	\$
<input type="checkbox"/> Bodypoint 1" Band Clamp Belt Mounting Kit	\$	IV Pole	\$
		<i>Not available with Transit Option or 8° Bend with Push Handle Back Post. Left side mount not available with O₂ Holder.</i>	
		<input type="checkbox"/> Right Side	<input type="checkbox"/> Left Side

Frame Color

Standard Colors			
Gloss Finish		Matte Finish	
<input type="checkbox"/> Black	N/C	<input type="checkbox"/> Army Green	N/C
<input type="checkbox"/> Candy Blue	N/C	<input type="checkbox"/> Charcoal Gray	N/C
<input type="checkbox"/> Candy Purple	N/C	<input type="checkbox"/> Iridescent Burgundy	N/C
<input type="checkbox"/> Candy Red	N/C	<input type="checkbox"/> Iridescent Navy	N/C
<input type="checkbox"/> Copper	N/C	<input type="checkbox"/> Shadow Black	N/C
<input type="checkbox"/> Silver	N/C		
Optional Colors			
Extended lead time			
Gloss Finish		Matte Finish	
<input type="checkbox"/> Blue Skies	\$	<input type="checkbox"/> Emerald Green	\$
<input type="checkbox"/> Bubblegum Pink	\$	<input type="checkbox"/> Retro Red	\$
<input type="checkbox"/> Granny Smith	\$	<input type="checkbox"/> Royal Blue	\$
<input type="checkbox"/> Metallic Orange	\$	<input type="checkbox"/> Smooth White	\$
<input type="checkbox"/> Teal	\$		
<input type="checkbox"/> Yellow	\$		

Cushions and Backrests



Cushion:

Backrest: