

# matrx™ pb

## POSTURE BACK

MaTRx-PB Posture Back is contoured to provide midline support and positioning. The extensive range of sizes provides exceptional comfort and support for a wide range of user needs.

### OUTER COVER:

- Moisture resistant and breathable
- Available in black or Kid\*ab\*ra print

### SACRAL/LUMBAR SUPPORT:

- Includes removable positioning pad

### BACK SHELL:


- Contoured shell (standard/deep) provides exceptional midline support and positioning
- Durable, lightweight ABS back shell includes built-in plate for headrest attachment

### FOAM:

- Molded foam wraps around shell for maximum comfort and protection

### HARDWARE:

- "Easy Set" hardware is easily installed and adjusted in minutes — adjustments can be made while user is seated in the wheelchair
- Quick release for easy portability
- Infinite adjustment of depth (4"), angle (40°) and rotation
- Fits 3/4", 7/8" and 1" diameter back canes
- Adjusts to fit chairs 1" narrower to 1" wider




Kid\*ab\*ra fabric optional



### 5" CONTOUR DEPTH-STANDARD



### 7" CONTOUR DEPTH-DEEP



### STANDARD POSTURE BACK

width	14"	16"	17"	18"	20"
height					
12"	●	●	●	●	●
16"	●	●	●	●	●
20"		●	●	●	●

Weight capacity 300 lbs.

### HEAVY DUTY POSTURE BACK

width	20"	22"	24"
height			
16"	●	●	●
20"	●	●	●

Weight capacity 500 lbs.

### EXTRA WIDE POSTURE BACK

width	24"	25"	26"	27"	28"	29"	30"
height							
16"	●	●	●	●	●	●	●
20"	●	●	●	●	●	●	●

Weight capacity 500 lbs.

### ➤ CANE REPLACEMENTS



- Replaces existing back canes which are less than 1" diameter for use with any MaTRx back

### ➤ FOAM-IN-PLACE KIT



- Individually molded insert produced in minutes with user seated in back
- Optional high stretch cover available

### ➤ EXTENDED MOUNTING PINS



- Extended mounting pins available for MaTRx-PB & Elite. Not available on Heavy Duty PB backs
- Allows the back to fit chairs up to 2" wider

### ➤ VISCO FOAM UPGRADE



- Laminate of 1/2" high resilient foam and 5/8" Visco Foam
- Available on MaTRx PB and Elite backs
- Molded foam is not available on MaTRx PB with this option

### ➤ HEADREST ADAPTER



- Allows the headrest to be mounted directly onto MaTRx Elite backs
- Available in 3 sizes