

## NAVIGAIT™

**NOW WITH FOOT PIECE TO BE USED ALSO INDOORS WITH NO SHOES OR SANDALS!**

NAVIGAIT™ is an alternative for users with a mild foot drop. The amount of dorsiflexion is easily adjusted with an elastic coil to accommodate the user's specific needs.

- Developed to help people with mild foot drop to improve dorsiflexion.
- The padded plastic shell sits comfortably on the front of the leg.
- Padding is washable.
- The amount of dorsiflexion can be adjusted.
- Only two sizes.

### **Intended Use**

NAVIGAIT™ is intended to support the foot when the ability to actively dorsiflex is reduced due to conditions such as Stroke, MS, Post Polio Syndrome.

### **Contraindications**

NAVIGAIT™ should not be used by patients who present spasticity or severe oedema. When used in combination with 4-FOOT, the fitter should pay extra attention to patients with diabetes.

### **Other**

NAVIGAIT™ does not address foot problems or deformities. These need to be treated with special shoes or insoles when needed.

For other indications we recommend the other products within the ToeOFF® Family.




As a new accessory to NAVIGAIT we now present 4-FOOT, a sock that makes it possible to use NAVIGAIT without shoes or shoes that do not have an ancor point (sandals). 4-FOOT comes in three sizes. Measure the circumference around the Metatarsal Heads for correct size.

4-FOOT has a plastic insert for increased stabilisation and "stay-up" silicon on the inside to keep it in place.



Item No.	Description	Size	Anterior section height
28410 0011	NAVIGAIT™	Small/Medium	27,5 cm
28410 0013	NAVIGAIT™	Large/XLarge	31,5 cm
28412 0011	Extra strap w/coil	Small/Medium	
28412 0013	Extra strap w/coil	Large/XLarge	
28413 0011	Extra strong strap w/coil	Small/Medium	
28413 0013	Extra strong strap w/coil	Large/XLarge	
			<b>Circumference</b>
28415 0011	4-FOOT Sock	Small	21 - 23 cm
28415 0012	4-FOOT Sock	Medium	23 - 26 cm
28415 0013	4-FOOT Sock	Large	26 - 29 cm



patent pending