

## 4 Wheel Folding Heavy-Duty Rollator

### Description

4 Wheel Folding Heavy-Duty Rollator

### Product Code

S20768



### Features and Benefits

- Lightweight with adjustable handle height to suit different users
- Extra width and strength for larger users
- Backrest and comfortable seat to allow use as a chair
- Simple folding mechanism for easy storage and transport
- Large diameter wheels with high-grip tyres and lockable brakes for safety
- Under-seat storage basket for convenience

## 4 Wheel Folding Heavy Duty Rollator

### Product Description

A lightweight but extra-wide and extra-strength rollator to help users who have difficulty or lack confidence when walking. The extra width and strength make this rollator suitable for larger users. Users can easily adjust the handle height to suit their needs. Brake levers on each handle control the brakes on the two rear wheels. They also include a locking mechanism for added safety. The large diameter high-grip wheels allow the rollator to roll along most surfaces while keeping control. A cushioned backrest and padded seat allow use of the rollator as a chair. The simple folding mechanism allows the user to easily collapse and store the rollator, useful when loading into a car boot or storage cupboard.

The rollator includes a removable under-seat storage basket to hold smaller items, keeping hands free for safety.

### Product Technical Details

- Painted metal frame with plastic parts.
- Lockable brakes to rear wheels.
- Removable storage basket.
- Maximum User Weight: 227 kg (35.5 st).
- Product weight: 9.4 kg.

