

REBOUND®  
POST-OP ELBOW  
Lightweight and simple to use




# REBOUND® POST-OP ELBOW

The Rebound Post-Op Elbow brace by Össur combines intuitive features for fitters and users with a design that aims to improve the post-op experience. Combining a drop-lock hinge and four telescoping paddles, individual post-op protocols are easily catered to. The Rebound Post-Op Elbow, lightest of its kind, promotes user comfort with straps that stay in place, anatomical left/right fit with moldable stability cuffs, and blue patient touch-points for easy education.

## Indications

Indications requiring elbow immobilisation and/or range-of-motion restriction. This may include:

- Stable fractures of the elbow, distal humerus, proximal radius or ulna
- Tendon and ligament injuries or repairs
- Collateral ligament reconstructions
- Elbow hyperextension
- Triceps tendon repair
- Chronic elbow injuries
- Tennis elbow release



### Simple and Intuitive

Easy education for users with intuitive patient touch-points, including Blue is You™ buckles and numbered strapping system. Telescoping paddles on the brace also enable individualised length adjustments to keep incision sites clear.



### Absolute Security

Drop-Lock hinge enables secure range-of-motion with bump-resistant pull adjustment stops and optional ROM lock-outs.



Light Weight

### Comfortable and Reliable

Designed to be the most lightweight brace on the market with strap guides that maintain strap location for added security and comfort.



## Ordering Information

Side	Part#
Left	B-8391680
Right	B-8391580

## Accessories

Description	Part#
Elbow Wrist Control Kit	B-8390003
Elbow ROM Locks Kit	B-8390005

## Replacement Part Kits

Description	Part#
Elbow Liner Kit	B-8390002
Elbow Sling Kit	B-8390004
Strap Tab Kit	B-8390006
Adhesive Hook Sticker Kit	B-8390007
Anti-Migration Pad Kit	B-8390008
Elbow Strap Kit	B-8390009



**ALLIED  
MEDICAL**  
LIMITED

**ÖSSUR**  
LIFE WITHOUT LIMITATIONS