



Walk Jog or Run – Have a Lot of Fun

Whether you search recreational leisure, physical challenge or just equipment that get's you from A to B, the PETRA RaceRunning bike does the trick!

Petra supports your body and sets your power free! On the track, in the forest or up and down the street. It is all about the joy of getting out! It is fun and stimulating regardless if you are a child, teenager or senior.

Your choice – Your pleasure.

Get on a PETRA and GO!

It is a fun and low impact activity that can be enjoyed with family, friends, non- or disabled peers.

It will open a whole new perspective on your abilities and opportunities well as a wide range of social interaction and learning experiences.

Same bike – same user – multiple purposes

- Basic mobility
- Play time
- Family rides
- Sport
- Physical education

We offer many alternative saddles and supportive parts to provide adequate stability and mobility on the PETRA.

With the spring assisted steering also people having major problems with spastisity balance and coordination can use PETRA.

It is easy to get from eg. your powerchair to PETRA because of the drop link saddle section.





PETRA Features

- Sizes from 3-4 years to adult.
- If the user walk "tip toe" – then add about 5 cm to leg length to find relevant saddle height.
- All sizes have saddle and chest support on multi adjustable fittings and posts.
- Each size can be adjusted app. 16 cm in saddle height and longer posts can be supplied.
- Adjustable handlebar with brake, two parking brakes on the frame and reflectors in wheels, rear and front.
- Reliable handling by assistive steering spring system.
- Mini PETRA, children age 3-6 has no drop link, usual equipped with junior bike saddle and PUR semi circle body support.
- Size X-small and up all have drop link of the saddle section.
- Durable & stabile design.
- Easy transfer on and off by use of drop link of saddle section.
- Tool bag and manual.


Measurements	0. Mini	1. X-Small	2. Small	3. Medium	4. Large	5. X-large
User max. height	110 cm	140 cm	140 cm	160 cm	185 cm	193 cm
User max. weight	30 kg	50 kg	50 kg	65 kg	85 kg	100 kg
Saddle height	38 - 53 cm	45-60 cm	50-66 cm	60 - 75 cm	70 - 85 cm	80 - 95 cm
Body support	50 - 75 cm	56-82 cm	63- 88 cm	75 - 102 cm	85 - 111 cm	95 - 125 cm
Handlebar	50 - 60 cm	60-70	65-75 cm	70 - 85 cm	80 - 95 cm	100 - 120 cm
Total length	118 cm	128 cm	144 cm	156 cm	168 cm	184 cm
Total width	71 cm	73 cm	75 cm	78 cm	82 cm	85 cm
Weight	10 kg	14 kg	14 kg	15 kg	17 kg	18 kg
Front wheels	17"	17"	20"	20"	20"	20"
Rear wheels	20"	20"	700cc	700cc	700cc	700cc



Belt and chest plate



Hand cuffs



2 PUR supports



Moon Plate behind



Droplink