

Posture Back Series: Sizing Chart

MaTRx Posture Back - contour depth 5"

Code	Width	Height	Fits wheelchair widths
MPB1412	14"	12" (short)	13" - 15"
MPB1416	14"	16" (regular)	13" - 15"
MPB1612	16"	12" (short)	15" - 17"
MPB1616	16"	16" (regular)	15" - 17"
MPB1620	16"	20" (tall)	15" - 17"
MPB1812	18"	12" (short)	17" - 19"
MPB1816	18"	16" (regular)	17" - 19"
MPB1820	18"	20" (tall)	17" - 19"
MPB2012	20"	12" (short)	19" - 21"
MPB2016	20"	16" (regular)	19" - 21"
MPB2020	20"	20" (tall)	19" - 21"

Deep Contour Version - contour depth of 7"

Code	Width	Height	Fits wheelchair widths
PBD1616	16"	16"	15" - 17"
PBD1620	16"	20"	15" - 17"
PBD1816	18"	16"	17" - 19"
PBD1820	18"	20"	17" - 19"
PBD 2016	20"	16"	19" - 21"
PBD2020	20"	20"	19" - 21"

Heavy Duty Version - contour depth of 5"

Code	Width	Height	Fits wheelchair widths
HDB2016	20"	16"	19" - 21"
HDB2020	20"	20"	19" - 21"
HDB2216	22"	16"	21" - 23"
HDB2416	24"	16"	23" - 25"

All Heavy Duty backs come standard with heavy duty hardware.