

Powerchair Script



Customer Details

Client _____
 Assessor _____
 Organisation _____

Client Weight Under 100kg Over 100kg* *Please specify _____

Base Type (Select PWC)	
Edge HD - Mid wheel, 204kg max user weight.	<input type="checkbox"/>
J623 2.0 - Mid wheel, 136kg max user weight.	<input type="checkbox"/>
Q1450 - Front wheel, 272kg max user weight.	<input type="checkbox"/>
R44 - Rear wheel, 136kg max user weight.	<input type="checkbox"/>

Tyres	
Pneumatic	<input type="checkbox"/>
Solid - Standard for Z, HD, 1450	<input type="checkbox"/>
Knobby Tyre	<input type="checkbox"/>

Powered Seating	
Non-Tilt	<input type="checkbox"/>
Tilt *	<input type="checkbox"/>
Lift n Tilt (REAC) ***	<input type="checkbox"/>

Power Recline	
Yes	<input type="checkbox"/>
No	<input type="checkbox"/>

Charger Required	
Yes	<input type="checkbox"/>
No	<input type="checkbox"/>

Seat Pan Width	
12"	<input type="checkbox"/>
13"	<input type="checkbox"/>
14"	<input type="checkbox"/>
15"	<input type="checkbox"/>
16"	<input type="checkbox"/>
17"	<input type="checkbox"/>
18"	<input type="checkbox"/>
19"	<input type="checkbox"/>
20"	<input type="checkbox"/>
21"	<input type="checkbox"/>
22"	<input type="checkbox"/>
23"	<input type="checkbox"/>
24"	<input type="checkbox"/>
Other - Please specify	<input type="checkbox"/>

Seat Pan Depth	
12"	<input type="checkbox"/>
13"	<input type="checkbox"/>
14"	<input type="checkbox"/>
15"	<input type="checkbox"/>
16"	<input type="checkbox"/>
17"	<input type="checkbox"/>
18"	<input type="checkbox"/>
19"	<input type="checkbox"/>
20"	<input type="checkbox"/>
21"	<input type="checkbox"/>
22"	<input type="checkbox"/>
23"	<input type="checkbox"/>
24"	<input type="checkbox"/>
Other - Please specify	<input type="checkbox"/>

Seat to Floor	
17.5"	<input type="checkbox"/>
18.5"	<input type="checkbox"/>
19.5"	<input type="checkbox"/>
Other - Please specify	<input type="checkbox"/>

Stump Support	
Left	<input type="checkbox"/>
Right	<input type="checkbox"/>
Both	<input type="checkbox"/>
Not Applicable	<input type="checkbox"/>

Rigidiser Bar	
Yes	<input type="checkbox"/>
No	<input type="checkbox"/>

Controller Type	
VR2 - J623 2.0 only	<input type="checkbox"/>
Q-Logic 3EX - Standard Hand Controller	<input type="checkbox"/>
Q-Logic 3E - Simplified Hand Controller	<input type="checkbox"/>
Attendant Control	<input type="checkbox"/>
Other - Please specify	<input type="checkbox"/>

Controller Side	
RHS	<input type="checkbox"/>
LHS	<input type="checkbox"/>
Midline	<input type="checkbox"/>
Not Applicable	<input type="checkbox"/>

Controller Mount	
Fixed In-line	<input type="checkbox"/>
Swing Away	<input type="checkbox"/>
Not Applicable	<input type="checkbox"/>

Armrest Type	
2 Post Flip Back	
Flip-up Armrests - Not compatible with recline	
Quick Height Adj Drop-In	
Heavy Duty	
Paediatric	
Recline Armrests - Only compatible with recline	
None	

Armrest Length	
Desk	
Full	

Armrest Length	
Standard	
Gel	
Channel Arms	

Legrest Type	
Centre Mount	
Powered Centre Mount *	
70 Degree Swing Away	
Powered 70 Deg Swing Away **	
60-90 Degree Swing Away	
Manual Elevating	
Angle Adjustable Centre Mount	
Powered C/Mount with Articulating Platform (Not available on J623 2.0)	
None	

Heel Loops	
Yes	
No	

Footrest Type/Size	
Standard - Std Only	
6 X 8 - AA Only	
8 X 11 - AA Only	
5 x 8 - Centremount Only	
5 x 11 - Centremount Only	
Paediatric	
Not Required	

- If Centremount selected, select footplate option below.

Seat Pan to Footrest - Left	
14°	
15°	
16°	
17°	
18°	
19°	
20°	
Other - Please specify	
Not Required	

Seat Pan to Footrest - Right	
14°	
15°	
16°	
17°	
18°	
19°	
20°	
Other - Please specify	
Not Required	

Footplate Options for Centre Mount:	
No F/Plate selected	
Bell Shaped F/Plate Length Ext	
Tapered F/Plate Width Ext	
Tapered F/Plate Outside Edge	
Tapered F/Plate Width Ext / Outside Edge	

Hip Belt Type	
Standard	
2 Point	
4 Point	
Other - Specify	
Not Required	

Hip Belt Size	
Small	
Large	
Medium	
Other - Specify	
Not Required	

Cushion Type	
Not Required	
Varilite Stratus	
Varilite Evolution	
Ottobock Cloud	
Ride Java	
Ki Mobility Axiom P	
Ki Mobility Axiom SP Visco	
Ki Mobility Axiom SP Fluid	
Other - Please specify	

Back Rest	
Icon Mid	
Icon Tall	
Icon Deep	
Ride Java Regular	
Ride Java Tall	
Ride Java Tall Deep	
Ride Java Grande	
Ride Java Grande Deep	
Other - Please specify	
Not Required	

Headrest	
Stealth 10" Non-Flip Down Mount	
Stealth 10" Flip Down Mount	
Stealth 14" Non-Flip Down Mount	
Stealth 14" Flip Down Mount	
Other (please specify)	
Not Required	

Laterals	
Yes	
No	

Tray	
Yes	
No	

Laterals Notes

Tray Notes

Specialty Controls
<p>Refer to Alternative Drive Controls Script for options</p>

Accessories

Set-Up Instructions