

EDGE 3



QUANTUM®

#1 FOR REHAB POWER

The Edge 3 is a revolutionary Mid-wheel drive power chair that provides a smooth and comfortable ride. Smooth ride suspension (SRS) is a standard feature, with iLevel as an option. iLevel is available at speeds up to 7.2kph with up to 12" of power adjustable seat height.



Features:

- Upgraded SRS (Smooth Ride Suspension)
- 7.2kph travel speed with iLevel activated
- Bold new colours



ALLIED
MEDICAL
LIMITED

SPECIFICATIONS




MODEL	Q6 Edge 3
DRIVE WHEELS	14" Pneumatic 14" Solid (available)
CASTOR WHEELS	Front: 6" Solid Rear: 6" Solid
MAXIMUM SPEED	Up to 10kph Up to 7.5kph with iLevel activated
GROUND CLEARANCE	2.7" / 6.9cm (battery tray)
TURNING RADIUS	20.5" / 52.1cm
OVERALL LENGTH	35.5" / 81.9cm without front riggings
BASE WIDTH	24.4" / 62.0cm
SEATING SIZES	Synergy/TRU-Balance 3 Static Seat: W: 12-22", D:12-22" TRU-Balance 3 Power Tilt: W: 12-22", D:12-22" TRU-Balance 3 Tilt and iLevel: W: 12-22", D:12-22"
SEAT TO FLOOR HEIGHTS	Synergy/TRU-Balance 3 Static Seat: 16.5-18.5" TRU-Balance 3 Power Tilt: 17.5-19.5" TRU-Balance 3 Tilt and iLevel: 17.5-19.5"
MANUAL RECLINE	Yes
BATTERY SIZE	55 Ah
BATTERY WEIGHT	17.2kg
AVAILABLE ELECTRONICS	120A Q-Logic 3 EX 75A Q-Logic 2 NE 75A Q-Logic 2 NE+
BATTERY CHARGER	8A off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel
WEIGHT CAPACITY	136kg
BASE WEIGHT	63kg

BENEFITS OF SRS

Smooth Ride Suspension is a superior automotive-grade powerchair suspension developed after two years of research and consumer feedback.

With unprecedented comfort and ride quality, SRS delivers exceptional driving performance over multiple terrains.

When you want a smooth, comfortable ride, unmatched performance and advanced stability, there's only one choice: SRS.



- Three coil over gas charged shocks provide an ultra-smooth and quiet ride, with exceptional performance over multiple terrains.
- Experiencing high amounts of vibration in a power wheel chair has been linked to user fatigue and discomfort. SRS technology can provide a comfortable ride, with unmatched performance and stability.
- A smooth ride can reduce environmental triggers for clients with spasticity.
- Significant decrease in the transmission of surface inconsistency
- Can significantly reduce pain levels, jolting forces are directly beneath the user.
- For clients with balance challenges, SRS can minimise the displacement of the upper body, head and neck. Providing a greater sense of stability and reducing postural requirements.



BENEFITS OF iLEVEL

- To transfer from surfaces (beds, toilets, shower chairs, car/van seats, etc.) that are higher than the seat on motorized wheelchairs
- To reach high enough to access one or more of the following areas at home, school, work, or the community:
 - Dresser drawers; clothes rod, washer or dryer
 - Medicine cabinet; bathroom sink/faucet; mirror; shower head/faucet
 - Freezer/refrigerator; oven/stove (including light/exhaust fan switches); microwave, cupboards/shelves; counter; sink/faucet; disposal switch
 - Light switches; thermostat; fire alarm; wall phone; speaker; door eye hole/viewer; doorbell; elevator buttons
- Difficulty seeing where you are going, moving through crowds, or safely crossing the street from a low (standard) seat height in motorized wheelchairs
- Reducing neck pain and difficulty seeing or hearing those standing nearby when seated in power wheelchairs

Contact your Business Territory Manager today, to experience the ride!