

Seahorse Plus Toileting & Shower Chair

The Seahorse Plus enables comfortable and safe toileting and bathing for children aged 5+ to adult



Designed & developed by NRS Healthcare

The Seahorse Plus inspires confidence and makes toileting and bathing a more manageable experience for both the carer and user.

Suitable for ambulant and non-ambulant users who require postural support, as well as users with mild to severe cognitive impairment or reduced motor function.

Fully adjustable for user security and comfort the Seahorse Plus can quickly and easily be positioned for use over a toilet or bath, as well as in the shower. It can also be used with a suitable hoist. Please contact the Customer Service team for advice on compatible hoists.

Available in three sizes: small, medium and large, with a range of comfortable PU liners to accommodate children as they grow, from 5 years to adult.

A range of optional accessories are available, including; liners, harnesses and anterior supports.

Product features

Tilt-in-Space



Tilt-in-Space: The innovative gas-assisted tilting mechanism enables 40° of rotation. To recline the seat simply press the red handle buttons then guide the user from a sitting position to fully reclined for safe, comfortable showering.



100mm seat height adjustment range enables optimal set up to suit most toilets and carer heights.





Comfortable removable armrests that can be substituted for the optional flip-up anterior supports.



Height and depth adjustable headrest with a lateral adjustment feature, to allow for variations in user head position.



Height adjustable footrests with soft heel guard and foot straps that feature a new lockable, flip up design to make entry and exit easy.



Multi adjustable lap-strap with central buckle and tightening feature.



Removable pommel / splash-guard aid user security and positioning.



Removable commode pan.

Seahorse Plus Optional Accessories

A range of optional accessories have been designed around the user and carer to offer additional safety, security and maximum comfort.



Seahorse Plus Liners:

Easy to fit, these soft PU foam liners have been designed for comfort and hygiene. The liners can be supplied in different thicknesses, allowing the chair to accommodate the user as they grow.

Small Thick Liner: N72310

Small Thin Liner: N72322

Medium Thick Liner: N72334

Medium Thin Liner: N72346

Large Extra Thick Liner: N84981

Large Thick Liner: N72668

Large Thin Liner: N72679



Chest Harness:

Optional Poziform™ Chest Harness for users requiring more secure positioning when using the Seahorse Plus Toileting and Shower Chair.

Small: N72139

Medium: N72140

Large: N72723



Butterfly Harness:

Particularly suitable for women the optional Poziform™ Butterfly Harness offers more secure positioning when using the Seahorse Plus Toileting and Shower Chair.

Small: N72152

Medium: N72164

Large: N72759