


# Seahorse Plus Dimensions



Feature	Ref	Small Seahorse Plus N72012			Medium Seahorse Plus N72024			Large Seahorse Plus N72036													
		Thick Liner N72310	Thin Liner N72322	No Liner	Thick Liner N72334	Thin Liner N72346	No Liner	Extra Thick Liner N84981	Thick Liner N72668	Thin Liner N72679	No Liner										
Seat Depth	SD	327mm	12¾"	355mm	14"	374mm	14¾"	425mm	16¾"	457mm	18"	470mm	18½"	455mm	18"	455mm	18"	455mm	18"	488mm	19¼"
Rear Width of Seat	SWR	258mm	10"	294mm	11½"	332mm	13"	303mm	12"	339mm	13½"	378mm	15"	375mm	14¾"	411mm	16"	447mm	17½"	477mm	19"
Front Width of Seat	SWF	294mm	11½"	334mm	13¼"	365mm	14¼"	353mm	14"	393mm	15½"	423mm	16¾"	420mm	16½"	460mm	18"	500mm	19¾"	530mm	20¾"
Height of Back	BH	391mm	15½"	391mm	15½"	414mm	16¼"	490mm	19¾"	490mm	19¾"	514mm	20¼"	549mm	21½"	549mm	21½"	549mm	21½"	572mm	22½"
Min Headrest height	HH	456mm	18"	456mm	18"	475mm	18½"	562mm	22"	562mm	22"	584mm	23"	620mm	24½"	620mm	24½"	620mm	24½"	641mm	25"
Max Headrest height	HH	631mm	25"	631mm	25"	650mm	25½"	737mm	29"	737mm	29"	759mm	30"	795mm	31½"	795mm	31½"	795mm	31½"	816mm	32"
Min height of Seat	SH	540mm	21½"	540mm	21½"	521mm	20½"	521mm	20½"	521mm	20½"	501mm	19½"	521mm	20½"	521mm	20½"	521mm	20½"	501mm	19½"
Max height of Seat	SH	648mm	25½"	648mm	25½"	629mm	24¾"	625mm	24½"	625mm	24½"	609mm	24"	629mm	24¾"	629mm	24¾"	629mm	24¾"	609mm	24"
Min height of Footrests	FH	275mm	10¾"	275mm	10¾"	255mm	10"	325mm	12¾"	325mm	12¾"	305mm	12"	325mm	12¾"	325mm	12¾"	325mm	12¾"	305mm	12"
Max height of Footrests	FH	375mm	14¾"	375mm	14¾"	355mm	14"	425mm	16¾"	425mm	16¾"	405mm	16"	450mm	17¾"	450mm	17¾"	450mm	17¾"	430mm	17"
Total width	OW	532mm	21"	532mm	21"	532mm	21"	586mm	23"	586mm	23"	586mm	23"	692mm	27"	692mm	27"	692mm	27"	692mm	27"
Total depth (upright)	OD	807mm	31¾"	807mm	31¾"	807mm	31¾"	860mm	33¾"	860mm	33¾"	860mm	33¾"	877mm	34½"	877mm	34½"	877mm	34½"	877mm	34½"
Max turning circle (tilted)	TC	1209mm	47½"	1209mm	47½"	1209mm	47½"	1413mm	55¾"	1413mm	55¾"	1413mm	55¾"	1488mm	58½"	1488mm	58½"	1488mm	58½"	1488mm	58½"
Min height over toilet	TH	461mm	18"	461mm	18"	461mm	18"	442mm	17½"	442mm	17½"	442mm	17½"	442mm	17½"	442mm	17½"	442mm	17½"	442mm	17½"
Max height over toilet	TH	569mm	22½"	569mm	22½"	569mm	22½"	551mm	21¾"	551mm	21¾"	551mm	21¾"	550mm	21¾"	550mm	21¾"	550mm	21¾"	550mm	21¾"
Min width over toilet	TW	426mm	16¾"	426mm	16¾"	426mm	16¾"	478mm	18¾"	478mm	18¾"	478mm	18¾"	592mm	23¾"	592mm	23¾"	592mm	23¾"	592mm	23¾"
Armrest to Seat	AH	179mm	7"	179mm	7"	199mm	8"	214mm	8½"	214mm	8½"	234mm	9"	214mm	8½"	214mm	8½"	214mm	8½"	234mm	9"
Seat angle	SA	0°-40°		0°-40°		0°-40°		0°-40°		0°-40°		0°-40°		0°-40°		0°-40°		0°-40°		0°-40°	
Max User Weight		64kg	10st	64kg	10st	64kg	10st	102kg	16st	102kg	16st	102kg	16st	152kg	24st	152kg	24st	152kg	24st	152kg	24st
Weight (including liners)		21kg	46.3lbs	20.2kg	44.5lbs	18.7kg	41.2lbs	26.1kg	57.5lbs	24.7kg	54.5lbs	22.3kg	49.2lbs	29.0kg	63.9lbs	28.3kg	62.4lbs	27.7kg	61.1lbs	24.8kg	54.7lbs