


NEW DESIGN

# Seat Mounting Brackets

Redesigned top-mount and **new** side-mount brackets provide an easy and secure way to mount belts to active wheelchairs.

Bodypoint® Seat Mounting Brackets — already your best mounting solution for chairs with sling seat upholstery or a solid seat base — are now even more versatile.



**New!** Side Mounting Bracket (20° bend)  
Part No. HW200-20-2 (pair)



Top Mounting Bracket (90° bend)  
Part No. HW200-90-2 (pair)  
(shown with saddle spacers, not included)