


# Special Tomato CatchUp Chair™



The seat that “catches up” with your child’s growth.

- Promotes proper body position and support
- One size seat fits most kids up to 113Kg.
- No-tool seat & footplate height & depth adjustments
- Small & Large Chest Harnesses provided
- Easy to clean & disinfect
- 5" locking casters for easy maneuverability.



The Special Tomato CatchUp Chair's unique design allows the seat depth to adjust from 11" to 17" and the seat to floor height adjusts 13" to 29" all with NO TOOLS! The compact design is excellent for use in schools, therapy clinics, hospitals, and home use.



*Comes standard with seat & back pads, pelvic harness, small & large chest supports. (shown)*

*Shown with optional size 3 seat & back liners, armrests and tray with storage bar*

## CatchUp Chair Specifications

| ITEM   | ENGLISH          | METRIC         |
|--|------------------|----------------|
| User Weight Capacity                             | 250 lbs          | 91 kg          |
| Armrest Pad- SMALL                               | 10"              | 25 cm          |
| Seat Width (Between Armrests) SMALL              | 14"              | 36 cm          |
| Armrest Pad- LARGE                               | 14"              | 36 cm          |
| Seat Width (Between Armrests) LARGE              | 16"              | 39 cm          |
| Seat Width (No Armrests)                         | 16"              | 41 cm          |
| Seat Depth                                       | 11" to 17"       | 28 cm to 43 cm |
| Seat Back Height                                 | *14" to variable | 34 cm          |
| Backrest Extension Height (off the seat cushion) | 14" to 28"       | 36 cm to 71 cm |
| Winged Headrest Height Range                     | 10" to 28"       | 25 cm to 71 cm |
| Shoulder Strap Height                            | 14" to variable  | 33 cm          |
| Seat to Footrest (variable)                      | 4" to 18"        | 11 cm to 46 cm |
| Seat to Floor                                    | 13" to 29"       | 33 cm to 74 cm |
| Seat to Back Angle                               | 100°             |                |
| Seat Angle to Floor                              | 5°               |                |
| Chair Weight                                     | 75 lbs           | 34 kg          |
| Push Handle Height                               | 34"              | 86 cm          |
| Overall Width                                    | 29"              | 74 cm          |
| Overall Height ( Top of Seatback to Floor)       | 42"              | 107 cm         |
| Overall Length                                   | 37"              | 94 cm          |
| Caster Size                                      | 5"               | 13 cm          |

| CATCHUP CHAIR & ACCESSORIES NUMBERS | DESCRIPTION   |
|-------------------------------------|---|
| 60000000                            | CatchUp Chair (Includes Seat and Back Cushion + (Lg & Sm 8 Point Harness) |
| 60020109                            | CatchUp Tray Kit & Storage Bar SMALL                                      |
| 60030109                            | CatchUp Adjustable Armrests (Required for Tray) - Small                   |
| 60020209                            | CatchUp Tray Kit & Storage Bar LARGE                                      |
| 60030209                            | CatchUp Adjustable Armrests (Required for Tray) - Large                   |
| 60040000                            | CatchUp Backrest Extension  |
| 60050009                            | CatchUp Winged Headrest (Requires Backrest Extension)                     |



**A.**

**A.** Tray with Storage Bar  
**B.** Armrests  
**C.** Backrest Extension  
**D.** Winged Headrest




**B.**



**C.**



**D.**



*Will fit up to a size 3  
 Soft-Touch® Sitter  
 (available in Lilac,  
 Aqua, or Gray)*