



# uniLink<sup>TM</sup>

## Allied Medical Product Guide







**ALLIED  
MEDICAL**  
LIMITED





**Stealth Products**



# UniLink | Allied Medical Product Guide

| Code           | Product Description                             | Image   | Status              | Added Detail  | Old Version  |
|----------------|---|---|---------------------|---|--------------|
| Head Solutions |   |   |                     |   |              |
| SPIDHBT300-U   | UniLink iDrive Pro Series Head Array            |    | Available for trial | i-Drive Bluetooth Pro Series Head Support with 3-Proximity Sensors, Occipital Pad (E2330), Swing Away Hardware for Left and Right Side (E1028), 2 Spot Pads, Egg Switch for reset mode change and the i-Drive Smart Box.<br><br>Add .PKMULTIQL3 for QL3 or DWR1111H053 FOR QL2. | SPIDHBT300-1 |
| SPIDHBT200-U   | UniLink iDrive Tri Head Array                   |   | Available for trial | i-Drive Bluetooth Tri-Array, Head Support (E2330), 3 proximity sensors, non-proportional, with egg switch for reset mode change and i-Drive Smart Box.  | SPIDHBT200-1 |
| SPIDHBT310-U   | UniLink iDrive Pro Series Head Array Sip & Puff |  | Special Order Item  | i-Drive Bluetooth Pro Series Head Support with Sip & Puff, includes: 3 Proximity Sensors, non-proportional, with egg switch for reset mode change with i-Drive Smart Box.   | SPIDHBT310-1 |
| SPIDHBT450-U   | UniLink iDrive Junior Curved Head Array         |  | Special Order Item  | i-Drive Bluetooth Junior (Standard Profile) Curved Head Array.  | SPIDHBT450-1 |








| Code        | Product Description   | Image   | Status              | Added Detail   | Old Version |
|-------------|---|---|---------------------|--|-------------|
| SPCP10CCFPM | UniLink Comfort Plus 10in Headrest                              |    | Available for trial | 10in (25.40cm) 17mm Ball Metric Contour Headpad - Actual dimensions 11 7/32in (28.51cm) wide x 5 3/4in (14.60cm) tall.<br>Use mounting SPUHB4D33200. | SPCP2510    |
| SPCP14CCFPM | UniLink Comfort Plus 14in Headrest                              |    | Available for trial | 14in (35.56cm) 17mm Ball Metric Contour Headpad - Actual dimensions 14 5/32in (35.94cm) wide x 5 3/4in (14.60cm) tall.<br>Use mounting SPUHB4D33200. | SPCP150     |
| SPACP150M   | UniLink Adjustable Comfort Plus 14in                            |   | Available for trial | 14in (35.56cm) Metric Headpad.<br>Use mounting SPUHB2D33200.   | SPACP150    |
| SPCP19CCFPM | UniLink Comfort Plus 19in Cool Core Flat Plate, Metric, E0955   |  | Available for trial | 19in (48.26cm) 17mm Ball Metric Contour Headpad - Actual dimensions 19 1/2in (49.53cm) wide x 5 1/8in (13.01cm) tall.                                |             |
| SPCP06CCFPM | UniLink Comfort Plus Nino 6in Headrest Flat Plate Metric, E0955 |  | Available for trial | Paediatric Headrest.<br>6in (15.24cm) 17mm Ball Metric Contour Headpad - Actual dimensions 6 1/4in (15.87cm) wide x 4 1/4in (10.76cm) tall.          | SPCPN350    |










| Code   | Product Description  | Image   | Status             | Added Detail  | Old Version |
|--------|--|---|--------------------|---|-------------|
| SP500M | Stealth Ultra Large Occipital Pad 6 X 3.75in Metric, E0955         |    | Special Order Item | Large Occipital Pad (Curved) - Labeled 6x3 - Real Dimensions 7 1/4in x 4 1/4in (18.42cm x 10.80cm) - Contour Depth 1in (2.54cm) - Metric.                 | SP500       |
| SP510M | Stealth Ultra Small Occipital Pad 4.5 X 2in Metric, K0108          |    | Special Order Item | Small Occipital Pad (Flat) - Labeled 4.5x2 - Real Dimensions 5 in x 2in (12.70cm x 5.08cm) - Metric.  | SP510       |
| SP600M | Stealth Ultra Suboccipital Pad 8.5 X 2.75in Metric, K0108          |  | Special Order Item | Standard Suboccipital Pad - Labeled 8 1/2 x 2 3/4 - Real Dimensions 7in x 3in - (17.78cm x 7.62cm) - Metric.<br>Maximum Neck Width 5in (12.70cm).         | SP600       |
| SP610M | Stealth Ultra Small Suboccipital Pad 6.375 X 2.25in Metric, K0108  |  | Special Order Item | Small Suboccipital Pad - Labeled 6 3/8 x 2 1/4 - Real Dimensions 6 1/4 x 2 1/2in - (15.88cm x 5.72cm) - Metric.<br>Maximum Neck Width 4 1/2in (11.43cm).  | SP610       |
| SP620M | Stealth Ultra X Large Suboccipital Pad 10.5 X 2.75in Metric, K0108 |  | Special Order Item | Extra Large Suboccipital Pad - Labeled 10 1/2 x 2 3/4 - Real Dimensions 9 1/2 x 2 3/4 - (24.13cm x 6.99cm) - Metric.<br>Maximum Neck Width 6in (15.24cm). | SP620       |



| Code           | Product Description  | Image   | Status              | Added Detail  | Old Version   |
|----------------|--|---|---------------------|---|---------------|
| SPUHB2D33200   | UniLink Headrest Hardware 17mm Ball - Vertical Bar Length 7in  |    | Available for trial | 7in (17.78cm) Height Adjustable, 1/2in (1.27cm) Square Area Vertical Bar.<br>Use for power recline.   | SPTWB485-17mm |
| SPUHB2D232FD   | UniLink Headrest Flip Down Mount 1/2in Sq 7.5in Vb             |    | Available for trial | Flip Down Hardware 7 1/2in (19.05cm) Height Adjustment, 1/2in (1.27cm) Square Area Vertical Bar.<br>Use for power recline.                      |               |
| SPUHB4D33200   | UniLink Headrest Hardware 17mm Ball - Vertical Bar Length 12in |  | Available for trial | 1/2in (1.27cm) Square 12in (30.48cm) of Height Adjustment Vertical Bar.   | SPTWB480-1    |
| PUHB4D33200-QC | UniLink Ultra 12in Sq - Vertical Bar Length 12in               |  | Special Order Item  | 1/2in (1.27cm) Square 12in (30.48cm) Height Adjustment Vertical Bar.<br>Does not include UniLink Ultra Link Set (SPUHBSU23200) sold separately. | SPTWB480-QFDM |
| SPUHBSU23200   | UniLink Headrest Stealth Ultra                                 |  | Special Order Item  | Set for UniLink QCR Head Solutions Hardware.  | SPSU980       |



| Code            | Product Description                      | Image   | Status              | Added Detail  | Old Version |
|-----------------|--|---|---------------------|---|-------------|
| Trunk Solutions |  |   |                     |   |             |
| SPULP535050320  | UniLink 3X5 Lateral Pad Neoprene         |    | Available for trial | 3in X 5in (7.62cm x 12.70cm) Lateral Pad, Neoprene, metric 1in (2.54cm) spacing.      | SPTWBL-3X5  |
| SPULP546050320  | UniLink 4X6 Lateral Pad Neoprene         |    | Available for trial | 4in X 6in (10.16cm x 15.24cm) Lateral Pad, Neoprene, metric 1in (2.54cm) spacing.     | SPTWBL-4X6  |
| SPULP535057120  | UniLink 3X5 Gel Lateral Pad Neoprene     |  | Available for trial | 3in X 5in (7.62cm x 12.70cm) Gel Lateral Pad, Neoprene, metric 1in (2.54cm) spacing.  |             |
| SPULP546057120  | UniLink 4X6 Gel Lateral Pad Neoprene 1in |  | Available for trial | 4in X 6in (10.16cm x 15.24cm) Gel Lateral Pad, Neoprene, metric 1in (2.54cm) spacing. |             |
| SPULP557057120  | UniLink 5X7 Gel Lateral Pad Neoprene 1in |  | Available for trial | 5in X 7in (10.16cm x 15.24cm) Gel Lateral Pad, Neoprene, metric 1in (2.54cm) spacing. |             |










| Code         | Product Description                    | Image   | Status              | Added Detail | Old Version |
|--------------|--|---|---------------------|--------------|-------------|
| SPULS25000-L | Stealth Short Swing Away Lateral - LHS |    | Available for trial |              | SPTWBNL-L   |
| SPULS25000-R | Stealth Short Swing Away Lateral - RHS |   | Available for trial |              | SPTWBNL-R   |
| SPULS25100-L | UniLink Swing Away Lateral LHS         |  | Available for trial |              | SPTWBL-L    |
| SPULS25100-R | UniLink Swing Away Lateral RHS         |  | Available for trial |              | SPTWBL-R    |









| Code        | Product Description                 | Image   | Status              | Added Detail   | Old Version |
|-------------|-------------------------------------|---|---------------------|----------------|-------------|
| SPUIMCNI780 | UniLink 7/8in Cane Clamp            |    | Available for trial | 7/8in (2.22cm) |             |
| SPUIMCNI100 | UniLink 1in Cane Clamp              |    | Available for trial | 1in (2.54)     |             |
| SPUIMPCI100 | UniLink Pelvic Thigh 1in Cane Mount |  | Available for trial | 1in (2.54cm)   |             |
| SPUIMCRI100 | UniLink 1in Crescent Clamp          |  | Special Order Item  | 1in (2.54cm)   |             |
| SPUIMCRI780 | UniLink 7/8in Crescent Clamp        |  | Special Order Item  | 7/8in (2.22cm) |             |





| Code            | Product Description                                    | Image   | Status              | Added Detail   | Old Version |
|-----------------|--|---|---------------------|--|-------------|
| SPUIMBKI090     | UniLink 90 Degree Bracket Lateral Mount                |    | Available for trial | 90° Lateral Mount Bracket for planar back surfaces as Stealth's TFB I-Back, Diego or TFB Flex Back.          |             |
| SPUIMBKI125     | UniLink 125 Degree Bracket Lateral Mount (TB3/TB4)     |    | Available for trial | 125° Lateral Mount Bracket for curves backs like Quantum Rehab TRU-Balance 3 or 4                            |             |
| SPUIMBKI135     | UniLink 135 Degree Bracket Lateral Mount (ADI, Q-Back) |   | Special Order Item  | 135° Lateral Mount Bracket for contoured backs like Stealth's ADI Backs, ADI Q-Back 2 and TFB Contour backs. |             |
| *SPUIPFSI100KIT | UniLink Lateral Full Surface Contact and 17mm Ball Kit |  | Available for trial | Kit code includes SPUIPFSI100 and SPBLL03. Both parts are needed for Full Surface Contact.                   | SPLTHW511   |
| SPUIMASI090     | UniLink 90 Degree All Season Bracket                   |  | Available for trial |  |             |
| SPUIMTNTB4      | UniLink Mount Track Nut for Quantum Powerchairs        |  | Available for trial | Track Mount for Quantum Rehab PWC Tracks.  |             |
| SPUIMBKITB4     | UniLink TB4 Lateral Mounting Kit                       |  | Special Order Item  |  |             |










| Code                        | Product Description                                | Image   | Status              | Added Detail  | Old Version     |
|-----------------------------|--|---|---------------------|---|-----------------|
| Upper Extremities Solutions |  |   |                     |   |                 |
| SPUEP544050320              | UniLink 4X4 Elbow Pad Neoprene                     |    | Available for trial | 4in (10.16cm) x 4in (10.16cm) Elbow Pad with Neoprene cover and Metric Threads.           | SPARMSES-4X4    |
| SPUEP544057120              | UniLink 4X4 Gel Elbow Pad Neoprene                 |    | Available for trial | 4in (10.16cm) x 4in (10.16cm) Gel Elbow Pad with Neoprene cover and Metric Threads.       | SPARMSES-4X4G   |
| SPUEP53005032R              | UniLink 3X10 Elbow Pad Neoprene Right              |   | Available for trial | Right 3in (7.62cm) x 10in (25.40cm) Elbow Pad with Neoprene cover and Metric threads.     | SPARMSES-3X10-R |
| SPUEP53005032L              | UniLink 3X10 Elbow Pad Neoprene Left               |  | Available for trial | Left 3in (7.62cm) x 10in (25.40cm) Elbow Pad with Neoprene cover and Metric threads.      | SPARMSES-3X10-L |
| SPUEP53005712L              | UniLink 3X10 Gel Elbow Pad, Neoprene, Metric Left  |  | Special Order Item  | Left 3in (7.62cm) x 10in (25.40cm) Elbow Gel Pad with Neoprene cover and Metric threads.  |                 |
| SPUEP53005712R              | UniLink 3X10 Gel Elbow Pad, Neoprene, Metric Right |  | Special Order Item  | Right 3in (7.62cm) x 10in (25.40cm) Elbow Gel Pad with Neoprene cover and Metric threads. |                 |



| Code           | Product Description                                    | Image   | Status              | Added Detail  | Old Version |
|----------------|--|---|---------------------|---|-------------|
| SPUETB4021FL-Q | UniLink Fixed Elbow Hardware LHS                       |    | Available for trial | Fixed Left Elbow Stop Hardware for Flip-Back Armrest        |             |
| SPUETB4021FR-Q | UniLink Fixed Elbow Hardware RHS                       |    | Available for trial | Fixed Right Elbow Stop Hardware for Flip-Back Armrest.      |             |
| SPUETB4213SL-Q | UniLink Swing Away Elbow Block For TB4 Left            |   | Available for trial | Left Swing Away Elbow Stop Hardware for Flip-Back Armrest.  |             |
| SPUETB4213SR-Q | UniLink Swing Away Elbow Block For TB4 Right           |  | Available for trial | Right Swing Away Elbow Stop Hardware for Flip-Back Armrest. |             |
| SPUESTDTRNFL   | UniLink 1in LHS Mount With Single Tracks Nuts Standard |  | Available for trial | 1in (2.54cm) Left Elbow Block Link Set Mount (Standard).    |             |
| SPUESTDTRNFR   | UniLink 1in RHS Mount With Single Tracks Nuts Standard |  | Available for trial | 1in (2.54cm) Right Elbow Block Link Set Mount (Standard).   |             |






| Code                       | Product Description                                | Image   | Status              | Added Detail  | Old Version   |
|----------------------------|--|---|---------------------|---|---------------|
| Pelvic and Thigh Solutions |  |   |                     |   |               |
| SPUTP546056100             | UniLink Pelvic Pad 4X6                             |    | Available for trial | 4in (10.16cm) x 6in (15.24cm) Pelvic/Thigh Track Pad, metric, 1in (2.54cm) Track Spacing. |               |
| SPUTP548056100             | UniLink Pelvic Pad 4X8                             |    | Available for trial | Pelvic/Thigh Track Pad 4in (10.16cm) x 8in (20.32cm) comes with Coolcore cover.           |               |
| SPUTP546050300             | Stealth 4X6, Neoprene Lat Pad, Metric, 1in Spacing |   | Available for trial | 4in (10.16cm) x 6in (15.24cm) Pelvic/Thigh Pad, metric, 1in (2.54cm) spacing.             | SPTWBADD-4X6  |
| SPUTP548050300             | Stealth 4X8 Pelvic/Thigh Pad, Metric, 1in Spacing  |  | Special Order Item  | 4in (10.16cm) x 8in (20.32cm) Pelvic/Thigh Pad, metric, 1in (2.54cm) spacing.             | SPTWBADD-4X8  |
| SPUTP540050300             | Stealth 4X10 Pelvic/Thigh Pad, Metric, 1in Spacing |  | Available for trial | 4in (10.16cm) x 10in (25.40cm) Pelvic/Thigh Pad, metric, 1in (2.54cm) spacing.            | SPTWBADD-4X10 |
| SPUTP54D050300             | Stealth 4X12 Pelvic/Thigh Pad, Metric, 1in Spacing |  | Special Order Item  | 4in (10.16cm) x 12in (30.48cm) Pelvic/Thigh Pad, metric, 1in (2.54cm) spacing.            | SPTWBADD-4X12 |
| SPUTP54D056100             | UniLink Pelvic Pad 4X12                            |  | Available for trial | Pelvic/Thigh Track Pad 4in (10.16cm) x 12in (30.48cm) comes with Coolcore cover.          |               |





| Code           | Product Description                         | Image   | Status              | Added Detail   | Old Version |
|----------------|---|---|---------------------|--|-------------|
| SPUTT7A1260000 | UniLink Thigh Guide Short Height Adjustable |    | Available for trial | Thigh Guide short 5in (12.70cm) height-adjustable and removable. |             |
| SPUTT6A1260000 | UniLink Thigh Guide Tall Height Adjustable  |    | Available for trial | Thigh Guide tall 7in (17.78cm) height-adjustable and removable.  |             |
| SPUTT6A116O2-L | UniLink Long Multiangle Thigh Guide LHS     |   | Available for trial | Long 7in (17.78cm) Multi-Angle Thigh Guide Left Side.            |             |
| SPUTT6A116O2-R | UniLink Long Multiangle Thigh Guide RHS     |  | Available for trial | Long 7in (17.78cm) Multi-Angle Thigh Guide Right Side.           |             |
| SPUTT23119O2-L | UniLink Fixed Adj Thigh Guide LHS           |  | Special Order Item  | Left Fixed One Size and Adjustable Thigh Guide.                  |             |
| SPUT800000-Q   | UniLink Thigh Guide Receiver Quantum Only   |  | Available for trial | Hip Guide Receiver for Quantum PWC's only.                       |             |



| Code            | Product Description                      | Image  | Status              | Added Detail  | Old Version  |
|-----------------|--|--|---------------------|---|--------------|
| Other           |  |  |                     |   |              |
| SPBLL03         | UniLink 17mm Ball Kit With Steel Inserts |   | Available for trial | This is the part required for the FSC (full surface contact). If this part also fits the laterals, it would negate the need for FSC lateral mounts. |              |
| SPIBLL17M       | UniLink Additional 17mm Ball             |   | Special Order Item  |   | SPTWB-17MMBK |
| SPUNILINK-KIT-1 | UniLink Kit in Backpack                  |  | Special Order Item  | Includes a tackle box, preconfigured head, trunk, elbow and pelvic/thigh assemblies with pads.  |              |