ALLIED MEDICAL

PIPER Comfort

Hybrid semi adjustable standard rehab-buggy



CE

- Robust modern buggy
- Versatile buggy for daily activities
- 2 sizes
- Unique backrest recline
- Lightweight aluminium frame
- Comfortable&durable
- Super compact foldable design
- Load capacity up to 60kg


patron

PIPER

PIPER Comfort is the new state of the art comfort rehab-buggy. This buggy is determined for short stay usage of children with special needs. The Piper Comfort features a unique backrest recline system with the adjustable calf rest allowing the user almost lie flat. With a range of postural positioning options, Piper is both supportive and comfortable resulting in reduced fatigue and safety for your child. Piper with its fully foldable undercarriage, robust but lightweight and compact aluminium frame, 5-point harness in standard, large adjustable footrest make it a versatile buggy for daily activities. Piper is available in two adaptable sizes. The variability of this buggy and the ability to attach a wide range of accessories allows the user's varying needs to be safely catered for.

Features



adjustable calf rest

backrest recline adjustable



large Footrest height adjustable



canopy with extension



one piece foldable handle



seat inlay

footplate booster



seat width adjustable by FLEXI-pads set

The seat to footrest length can be adjusted and by using an optional additional seat size-reduction inlay and footplate booster the buggy can be adjusted to fit the child's changing size from a child of approximately two up to a maximum load of 50kg's.



front wheels with direction lock




super compact foldable design

sizes and patterns


sizes and patterns

2+ sizes






Denim blue



Teal green




Grey-black



Antique pink

indicative size guide




height of the user

with inserts	34	40
85 - 110 cm	95 - 130 cm	115 - 160 cm




maximum user weight

with inserts	34	40
50,0 kg	50,0 kg	60,0 kg



Olive green



Tabasco brown

colour range

colour range



Denim blue



Teal green



Grey black



Antique pink



Olive green



Tabasco brown

specification

Standard equipment:

- seat depth non-adjustable
- extended backrest height non-adjustable
- backrest recline adjustable
- footrest height adjustable, detachable
- seat width adjustable by FLEXI-pads set (optional eq.)
- 5-point immo.harness with crotch belt and shoulder padding
- 175/225 (175/265) mm solid wheels, front with swivel lock
- detachable wheels, forks
- one-piece foldable handle
- basket, transportation lock
- frame colours and cover patterns according to the current offer

The product had been successfully crash tested and fulfilled all the requirements of EN ISO 7176/19.
It enables to use the rehab-buggy not only for a common duty, but as well as for the transportation in specially adjusted motor vehicles, equipped with the anchoring points and equipment.

order codes

Products	
PIPER Comfort 34:	PIF34RP01
PIPER Comfort 40:	PIF40RP01

accessories



Canopy CLASSIC



Headrest adjustable



Abduction Block



Footplate booster

.... and much more accessories!

according to the actual list

	with Seat inlay / footrest raiser	34	40
Seat width [mm]	230/270	260/step/340	320/step/400
Seat depth [mm]	240	290	340
Backrest height [mm]	680	680	780
Seat to footplate [mm]	200-230	240-370	340-470
Seat angle	+25°	+25°	+25°
Backrest angle		90° - 125°	90° - 125°
Weight		14,5 kg	17,5kg
Maximum load		50,0 kg	60,0kg
Total length [mm]		1110	1310
Total width [mm]		590	635
Total height [mm]		1000	1060
Brake		foot oper.	foot oper.
Folded stroller length [mm]		1120	1220
Folded stroller width [mm]		330	440
Folded stroller height [mm]		280	450

notes